Endopuls 811

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Manual del usuario
Istruzioni per l'uso

Copyright:


Enraf-Nonius B.V.
P.O. Box 12080
Vareseweg 127
3004 GB Rotterdam
The Netherlands
Tel: +31 (0)10 – 20 30 600
Fax: +31 (0)10 – 20 30 699
info@enraf-nonius.nl
www.enraf-nonius.com

Part number: 1650751_40
November 9, 2012

Endopuls 811

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Manual del usuario
Istruzioni per l'uso

[1] LCD-display:

Kleurendisplay met touchscreen

[2] Controller links:

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1	Foreword	40
2	Product liability	40
3	Precautionary Instructions	40
4	Intended Use	41
5	Product description	41
6	Indications/Contraindications	41
	6.1 Indications:.....	41
	6.2 Contraindications:.....	42
	6.3 Precautions and Warnings:	42
7	Package Contents	43
8	Application Notes	44
	8.1 Shockwave	44
	8.2 Hand-piece	44
	8.3 Standby mode on device and hand-piece	44
	8.4 Applicators	44
	8.5 Changing applicators	45
	8.6 Foot switch.....	45
9	Operator controls and accessories	46
10	Installation	48
11	Basic Operation	49
	11.1 Home menu.....	49
	11.2 System Settings.....	50
	11.3 SD CARD	50
	11.4 Clinical Protocols	51
	11.5 Favorites.....	52
	11.6 Delete Favorites	54
	11.7 Manual Operation	55
12	Maintenance and Troubleshooting	57
	12.1 User Maintenance	57
	12.1.1 <i>Cleaning and disinfection of the apparatus</i>	57
	12.1.2 <i>Monitoring the hand-piece temperature</i>	57
	12.2 Troubleshooting.....	58
	12.2.1 <i>Failure or malfunction of the hand-piece</i>	58
	12.2.2 <i>Irregular delivery of shockwaves / overheating of hand-piece</i>	58
	12.2.3 <i>No response at main switch / display remains dark</i>	58
	12.3 End of life.....	58
13	Function Test	59
14	Error messages	59
15	Specifications	60
16	Safety and Performance standards	60
17	EMC details	61

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1 Foreword

This manual has been written for the owners and operators of the Endopuls 811. It contains general instructions on operation, precautionary practices, maintenance and parts information. In order to maximize the use, efficiency and lifespan of your unit, please read this manual thoroughly and become familiar with the controls as well as the accessories before operating the unit. Specifications put forth in this manual were in effect at the time of publication. However, owing to Enraf-Nonius her policy of continual improvement, changes to these specifications may be made at any time without obligation on the part of Enraf-Nonius.

2 Product liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortcomings of the product.

To the maximum extent permitted by applicable law, in no event will Enraf-Nonius or its suppliers or resellers be liable for any indirect, special, incidental or consequential damages arising from the use of or inability to use the product, including, without limitation, damages for loss of goodwill, work and productivity, computer failure or malfunction, or any and all other commercial damages or losses, even if advised of the possibility thereof, and regardless of the legal or equitable theory (contract, tort or otherwise) upon which the claim is based. In any case, Enraf-Nonius's entire liability under any provision of this agreement shall not exceed in the aggregate the sum of the fees paid for this product and fees for support of the product received by Enraf-Nonius under a separate support agreement (if any), with the exception of death or personal injury caused by the negligence of Enraf-Nonius to the extent applicable law prohibits the limitation of damages in such cases. Enraf-Nonius cannot be held liable for any consequence resulting from incorrect information provided by its personnel, or errors incorporated in this manual and / or other accompanying documentation (including commercial documentation) The opposing party (product's user or its representative) shall disclaim Enraf-Nonius from all claims arising from third parties, whatever nature or whatever relationship to the opposing party.

3 Precautionary Instructions

In this section general Warnings and Precautions are listed, that you should be aware of when using the Endopuls 811. See also chapter 6.3 for Warnings and Precautions that are application specific.

WARNING:

- Federal law (USA only) restricts this device to sale by, or on the order of, a physician or licensed practitioner. This device should be used only under the continued supervision of a physician or licensed practitioner.
- Make certain that the unit is electrically grounded by connecting only to a grounded electrical service receptacle conforming to the applicable national and local electrical codes.
- Do not operate the unit in an environment of short-wave or micro-wave diathermy.
- Care must be taken when operating this equipment around other equipment. Potential electromagnetic or other interference could occur to this or to the other equipment. Try to minimize this interference by not using other equipment in conjunction with it.
- This equipment is not suitable for use in the presence of flammable anesthetics mixture with air, oxygen, or nitrous oxide.
- This device should be kept out of the reach of children.

CAUTION:

- Read, understand and practice the precautionary and operating instructions. Know the limitations and hazards associated with using any electrical stimulation device. Observe the precautionary and operational decals placed on the unit.

- This unit should be operated in temperatures between 10 °C and 25 °C (50 °F and 77 °F), with a Relative Humidity ranging from 20%-80%.
- Do not expose the unit to direct sunlight, heat radiated from a heat radiator, excessive amounts of dust, moisture, vibrations and mechanical shocks.
- In the case of ingress of liquids, unplug the unit from the mains supply and have it checked by an authorized person (see the paragraph on maintenance)
- Before administering any treatment to a patient you should become acquainted with the operating procedures for each mode of treatment available, as well as the indications, contraindications, warnings and precautions. Consult other resources for additional information regarding the application of shockwave therapy.

4 *Intended Use*

The Endopuls 811 is a state of the art innovative shockwave treatment device. The radial, ballistic shockwave treatment is a procedure with a wide range of applications ranging from superficial orthopaedic problems to myofascial trigger point treatment.

The device is designed to only be used by or under the supervision of persons using the medical device in the course of their work and in the framework of a professional healthcare activity, who understand the benefits and limitations of shockwave therapy.
I.e. "professional users".

5 *Product description*

What does Endopuls 811 do?

Creation of shockwaves using an ergonomic hand-piece and the transmittal of the shockwaves via special applicators. The Endopuls 811 has a maximum penetration depth of about 35 mm in human tissue.

The technology:

An electromagnetic field is generated via a coil in the back of the hand-piece.
A projectile is accelerated as a result of the field; this strikes against the applicator at the front of the hand-piece and generates shockwaves, which spread out radially in the tissue.

The advantage:

The innovative technology allows a compact design with no need for a compressor.
The clear and modern colour display shows all relevant parameters for treatment and the modern touch operation ensures pleasure and motivation when providing treatment.
Individual programme start configuration and clear, simple menu navigation make operation of the device easy and comfortable for users.
Infinitely variable frequencies and various applicators allow treatment to be adapted to the particular condition of the patient.
The compact design saves room in the practice and is highly suited for home visits.

Note: The device should only be used by medical specialists (e.g., doctors, therapists and health paraprofessionals).

The Endopuls 811 has been constructed and designed solely for the treatment of superficial orthopaedic problems in humans and animals.

6 *Indications/Contraindications*

6.1 Indications:

- Radial and ulnar epicondylitis
- Calcific tendonitis of the shoulder / shoulder problems
- Status post muscular injuries
- Chronic patellar tendonitis

- Jumper's knee
- Achillodynia
- Plantar fasciitis
- Heel spurs
- Myofascial trigger point treatment e.g. neck
- Myofascial trigger point treatment e.g. back, muscular back pain
- Bursitis Trochanterica
- Periostitis / shin splints (status post strain)

6.2 Contraindications:

- vascular diseases present in or near the treatment area
- local infections in the treatment area
- around malignant or benign tumours
- directly on cartilage surfaces or near the small facet joints of the spinal column
- directly over implanted electronic devices such as pacemakers, analgesic pumps, etc.
- in areas, in which mechanical energy in the form of vibrations may lead to tissue damage such as metal implants after a fracture

In general we advise against treatment

- if blood clotting disorders are present or the patient is receiving treatment that results in a change in the blood clotting behaviour
- during pregnancy
- on patients with neurological diseases resulting in impairment of the vasomotor function in the treatment area
- over air-filled cavities such as treatment on the thoracic spine, etc.
- on children, particularly around the epiphyseal plates

Care is required for patients

- with impaired sensibility
- with severe autonomic disorders
- under the influence of drugs and/or alcohol as circulatory stresses and inadequate treatment responses cannot be excluded.

6.3 Precautions and Warnings:

- Users of the Endopuls 811 shockwave treatment device must be trained in how to use the system properly and have the appropriate skills.
- Any treatment instructions regarding treatment location, duration and strength require medical knowledge and should only be given by authorised doctors, therapists and health paraprofessionals. It is imperative that these instructions are followed.
- Treatment must always be carried out under medical supervision.
- The Endopuls 811 hand-piece is not designed for permanent use. After a treatment with max. 6000 shocks, a break of 15 min. becomes necessary.

Caution:

Patients who are concurrently receiving treatment involving a reduction and/or modification of blood clotting or prolongation of the blood clotting time (e.g. with acetylsalicylic acid) should consult their therapist about possibly stopping this treatment as these patients may be more prone to greater haemorrhaging and bruising when radial shockwaves are applied.

Shockwaves are strongly scattered in air pockets and create reflections that may have negative effects. You must therefore never perform any direct treatments over the lungs (intercostal spaces) or the gastrointestinal area.

It must not be used in wet areas. If it is used in wet areas, significant damage may occur, and patients and users may be endangered.

The instruments must only be operated with the mains cable provided. Protect the mains cable from any mechanical stress.

7 Package Contents

1650.900 Endopuls 811

Standard Accessories

1650800 Hand-piece, complete with a 15 mm applicator
1650808 Holder for hand-piece
1650801 25 mm applicator
1650802 15 mm applicator
1650803 6 mm applicator
1650804 Silicon covers, 10 pcs
1650805 Endopuls lotion, bottle 250 ml
1650806 Foot switch
1650807 Applicator measurement tool
1650809 Carrier case
1650751 Operating manual
xxxxxxx Mains cable

8 Application Notes

8.1 Shockwave

The Endopuls 811 operates with mechanical energy. The energy is transmitted to the patient via a hand-piece, which is usually held in one hand.

To do this, the hand-piece is placed on the area or point of treatment with the applicator held vertically. When the shockwave is activated, it is possible to work either steadily on a single spot or dynamically over an area.

It is advisable to use gel (included in the accessories) in order to reduce friction on the skin.

Due to the weight of the hand-piece it is normally not necessary to apply pressure to the treatment area / point.

The hand-piece is placed on the treatment area / point and is held loosely in position with one hand.

If required, additional pressure may be applied in the direction of the tissue, and the working angle can be varied.

CAUTION:

The use of water based gels or other lubricants are forbidden, oxidation could damage the hand-piece. We advise to use only the original lotion or only non-water based gel and lotions. The applicator must be covered with a silicone cap to protect it.

NOTE:

Despite high internal absorption as a result of the weight and design of the hand-piece, vibrations may cause strain to the user's hand.

- Recommended protective measure: Limit the duration of exposure
- The patient should be carefully monitored throughout the treatment.

8.2 Hand-piece

The hand-piece contains the shockwave generator, a fan to dissipate heat and the slot for the different applicators. It is connected to the control unit.

NOTE:

- The shockwave generator in the hand-piece is an expendable part and has to be replaced after a specific period of use, as its functionality decreases over time.
- Enraf-Nonius BV. guarantees unrestricted use of at least 2 million shocks per shockwave generator.
- Wear on the shockwave generator varies. Depending on performance and frequency, sometimes far more than 2 million shocks can be delivered.
- For more information on the need to replace the shockwave generator, see troubleshooting

CAUTION:

To work with the hand-piece on a patient, it is essential that one of the applicators is screwed tightly onto the hand-piece as far as it will go.

The cable should not be stretched beyond its maximum length and must be protected against pinching or any other mechanical damage.

To avoid heat accumulating in the hand-piece, it is essential to ensure that the hand holding it or anything else does not block the air vents at the top, and particularly, on the base of the hand-piece.

8.3 Standby mode on device and hand-piece

The fan in the hand-piece is activated by depressing the footswitch and stops automatically after reaching a certain temperature.

8.4 Applicators

There are 3 different applicators available for treatment.

8.5 Changing applicators


To change the various applicators, hold the hand-piece in one hand and unscrew the applicator from the hand-piece with the other hand (anticlockwise). Screw the required head tightly onto the hand-piece (clockwise), until the black outside ring of the applicator rests on the hand-piece (there should not be any thread visible).

NOTE:

- Applicators are expendable parts and must be replaced after a certain period of use.
- Minor / slight deformation or shortening of the rear impact dome does not affect functionality.
- In cases of greater deformation or stronger shortening of the rear impact dome the applicator must be replaced.

A test template is supplied with the device that enables the user to test if the wear limit has been reached (see diagram).


Applicator measurement

Air gap between applicator and the applicator measurement tool

Applicator is OK

Applicator measurement tool

Template makes contact or an air gap at tip

Applicator need to be changed


8.6 Foot switch

Place the foot switch so that it can be reached easily during treatment. The footswitch control unit is multi-directional so it is not necessary to align the footswitch exactly.

To avoid damage, please note that only slight pressure needs to be put on the switch. Use the front part of your foot, not the heel to operate the footswitch.

The switch does not have a locking device, which means that it only remains actuated as long as pressure is applied.

9 Operator controls and accessories


[1] LCD display:
Colour display with touch screen

[2] Controller left:
This controller can be used for the adjustment of the energy

[3] Controller right:
With this controller you can adjust the frequency and target.

[4] SD card slot:
This is the sleeve to insert the SD card.

[5] Mains Switch
To switch the unit on and off

[6] Connector for mains cable
To disconnect the device from the mains supply, unplug the power cord.

[7] Connection for the hand-piece
Hand-piece connection

[8] Connection for the foot switch

This is the connector to connect the foot switch

CAUTION:

- Connection of accessories other than the ones specified by the manufacturer can adversely affect the safety of the patient and correct functioning of the equipment, and is therefore not permitted.

[9] Type number/warning sticker

Provides information on the apparatus, such as type and serial number, as well as connection data such as mains voltage and maximum current consumption.

[10] Hand-piece:

Hand-piece with changeable applicators

[11] Hand-piece holder:

Can be positioned on the right or left side of the unit.

[12] Applicator 20mm

Applicator with large surface

[13] Applicator 15mm

Applicator with medium size surface

[14] Applicator 6mm


Applicator with small surface

[15] Silicone cap

Protects the applicator and is easy to clean

[16] Foot-switch:

To start the treatment. Foot-switch is placed on the floor

10 Installation

NOTE:

- Before starting up the system, remove the Endopuls 811 from its transport case. Do not operate the device while it is in the case.
- Ensure that the Endopuls 811 is placed on a stable surface
- Make sure that the main switch on the device is set to '0'.

Connecting the mains cable

- Insert the mains cable into socket [6] and connect it to a wall socket.

Connecting the hand-piece

- Plug the hand-piece into the appropriate socket [7] of the device and place it on the table

NOTE:

Ensure that an applicator is inserted into the hand-piece and that it is properly screwed in as far as it will go.

Connecting the footswitch

- Plug the footswitch into the appropriate socket [8] of the device and then place it on the floor.

Switching on the device

- Switch on the device using the main switch [5].


11 Basic Operation

Turn on the apparatus

- Switch on the device using the main switch [6].

11.1 Home menu


The Endopuls 811 is equipped with touch screen technology. Via the menu you have access to all functions. Select the desired function by pushing the button or icon.

NOTE:

Changes to the default settings can only be made from the start screen. Press button "System Settings" [A] to open the screen

11.2 System Settings

- Select your language via the ▲▼ icon Push ✓ to confirm and to return to the previous menu
- Adjust the brightness via the ▲▼ icon Push ✓ to confirm and to return to the previous menu
- Adjust the Volume via the ▲▼ icon Push ✓ to confirm and to return to the previous menu

Hand-piece Counter Status:

The counter status for the hand-piece that is currently connected, is shown in this display field.

Version info:

Press the version info button to open the window with information about the current software version of the device.

Restore defaults:

Press the restore default button to reset the factory default settings.

Touch Calibration:

Press the "Touch Calibration" button to open the screen to carry out the touch calibration.

This can be done to improve the touch input if it is not sufficiently accurate.

First press the + symbol in the top left corner. A + symbol will appear in the lower right corner.

Then precisely press the + symbol in the lower right corner.

Repeat the procedure to complete the touch calibration.

11.3 SD CARD

User-defined settings and the treatment recommendation list are saved on the SD card.

NOTE:

If the SD card is not inserted, the message 'SD card not found' appears when the Clinical Protocols and Favorites buttons are pressed.

Deactivate the message by pressing the button "OK" and continue.

11.4 Clinical Protocols


The Clinical Protocols are treatment recommendations. This helps you to select the treatment. The treatment can be selected from the treatment recommendation list. Push the Clinical Protocols button to access.

Select in the list the treatment recommendation by pushing that button or:


- Push icon Info [B] for more therapy information about the protocol or,
- Push the ▲▼ icon [C] to scroll in the list

Treatment information

- Push the ▲▼ icon [E] to select the next page.
- Push the accept button ✓ [F] to select the program, the treatment screen appears.

- The unit is ready to use.

NOTE:

It's always possible to change the parameters.


- Push the store button  to add this program in your favorites.


11.5 Favorites

Here you can store you favorite setting by name. The list is empty and you need to store some settings first.

Use manual settings or a treatment recommendation.


- Set the desired parameters and push the store button [F]

- Select a position in the Favorites list by pushing [A]
- Scroll in the list with the▲▼ icon [C]


- Insert a name via the keyboard


- Push the accept button [F] to store this program in the favorites list.

11.6 Delete Favorites


- Push the trash button [F]

- The button changes to a different color


- Select the position you want delete from the list [A]
- An pop-up appears


- Push the accept button [V] to delete the program
- Push the cancel button [X] to cancel the deleting procedure.

11.7 Manual Operation

Here you have immediate access to all parameters.

Push the button [A]“Manual Operation” in the home menu.

- Adjust the frequency by pushing the frequency read-out [G] and change the frequency with the right controller.
- Adjust the target by pushing the target read-out [H] and change with the right controller.
- Adjust the energy by pushing the energy read-out [I] and change with the left controller.
- Start the therapy via the footswitch.

Reset counter: push the RESET COUNTER button to reset to zero

Use the count up / countdown button to toggle between these two functions

Use the pulsed / continuous button to toggle between continuous and pulsed mode

NOTE:

If pulsed mode is selected you can select the 3 Burst modes.

12 Maintenance and Troubleshooting

12.1 User Maintenance

Separate servicing is not required for this product.

Before starting any maintenance or cleaning, the device must always be switched off at the main switch and the plug pulled out.

You should also check the applicator domes for any wear, as described in chapter 7.1.4.

WARNING:

When using lubricants, it is essential to pull the silicone cap over the applicator.

If you do not use the protective cap, the lubricant can get inside the applicator and hand-piece, which can lead to permanent soiling and malfunctioning.

NOTE:

In this case the warranty becomes void.

12.1.1 Cleaning and disinfection of the apparatus

Clean the device and hand-piece with soap lotions or cleaning agents that do not contain alcohol or solvents.

Conventional disinfecting products used for medical equipment are suitable. To disinfect the device, hand-piece and accessories we advise Dismozon®Pur. in a solution between 0,25 and 1,5%. (see at <http://www.bode-chemie.com> for the local supplier)

CAUTION:

- It is essential to ensure that no moisture gets into the system when cleaning.
- Do not spray the cleaning agent directly on the glass panel.
- Do not use cleaning agents that contain strong alkalis, lye, acid, detergents with fluoride or detergents with ammonia.

12.1.2 Monitoring the hand-piece temperature

Generating mechanical shockwave energy causes a considerable build-up of heat in the hand-piece.

To avoid shortening the lifetime of the hand-piece, a temperature switch has been integrated. This triggers an internal switch-off, if the temperature becomes too high, forcing the hand-piece to cool. The warning message appears at the screen.

As soon as the hand-piece has reached the correct operating temperature, automatically the pop-up will disappear and the unit is again ready for use.

If the temperature switch is activated, this is indicated by a message on the display and shocks can no longer be emitted.

12.2 Troubleshooting

12.2.1 Failure or malfunction of the hand-piece

- Check to ensure that the hand-piece plug is properly connected to the device. It must be fully engaged.
- Check the cable of the hand-piece for any mechanical damage

12.2.2 Irregular delivery of shockwaves / overheating of hand-piece

Possible cause 1: Wear of applicator. Difficult to move due to wear. Applicators are wear parts and should be replaced after a specific number of shocks.

Remedying cause 1:

- Removal of parts subject to abrasion:
- Remove the applicator from the hand-piece and clean the rear dome thoroughly. Then hold the hand-piece, without the applicator, with the opening downward and, at 2 or 5 Hz frequency, release a few shocks (maximum 10) at the lowest energy level. Then reinsert the applicator.
- If the error still occurs, the applicator has to be changed.

Possible cause 2: Wear of shockwave generator

The shockwave generator is an expendable part and should be replaced after 2 million shocks. Check the total number of shocks of the device in the configuration menu.

Remedying cause 2:

- If the total number of 2 million shocks has been reached or exceeded, the shockwave generator must be replaced.
- To replace the shockwave generator, contact your local distributor or a qualified customer engineer.

12.2.3 No response at main switch / display remains dark

- Make sure that the mains plug is properly inserted in the power outlet and the device connector is firmly plugged into the device port.
- Inspect the mains cable for damage.
- Check the power supply and the power plug.
- Above the mains input socket of the device, there are fine-wire fuses, which isolate the mains voltage in the event of an electrical problem. Open the flap and check the fuses.
- Replace any faulty fuses.

NOTE:

Only replace a fuse with one of exactly the same name or one that is equivalent. Before doing this, check the entire power supply for any possible faults.

If the error occurs again, it is essential to inform your local distributor.

12.3 End of life

The Endopuls 811 contains materials that can be recycled and/or are noxious to the environment. Specialized companies can dismantle the unit and sort out these materials. When you dispose of the unit, find out about local regulations concerning waste management.

13 Function Test

The Endopuls 811 will run a self-test that checks all internal components after it is switched on.

- An error message is shown in case of faults.
- In addition, a function test shall be made monthly or in case of doubt of the proper functioning of the device.

NOTE:

- Before performing the function test, check whether the hand-piece and the foot switch are connected correctly to the device.
- Check for proper mains connection.

Switch on the device.

Depress the footswitch briefly – the fan and shockwave generator will start immediately, whereby the shockwave generator has to operate at the frequency indicated on the display (5 Hz as default value).

NOTE:

In conclusion of the test, switch off the device at the main switch.

If a treatment is to be performed immediately afterwards, set the required treatment parameters and proceed as mentioned in Chapter 11.7

14 Error messages

Monitoring the hand-piece temperature (sww 12.1.2)

Generating mechanical shockwave energy causes a considerable build-up of heat in the hand-piece. To avoid shortening the lifetime of the hand-piece, a temperature switch has been integrated. This triggers an internal switch-off, if the temperature becomes too high, forcing the hand-piece to cool.


If the temperature switch is activated, this is indicated by a message on the display and shocks can no longer be emitted.

As soon as the hand-piece has reached the operating temperature, automatically the pop-up disappears and the unit is ready to use again.

No SD-Card found

If the SD card is not inserted, the message 'SD card not found' appears when the 'Favorites' and 'Memory' buttons are pressed.

Insert card and confirm with 'OK'.

15 Specifications

Endopuls 811	Treatment system for electromagnetic generation / application of of radial shockwaves in orthopaedics and physiotherapy.
Dimensions	L 350 mm / W 250 mm / H 140 mm
Weight	2.7 kg
Power supply	100–240 VAC / 50/60 Hz, 220 VAC / 60 Hz
Fuse	3,15 AT
Conformity	Protection class I / Application class BF
Frequency range	1 Hz – 22 Hz, can be adjusted in 1 Hz steps
3 burst modes	4, 8, 12 shots
Shock energy levels	60 – 180 in steps of 10 mJ and 185 mJ (at the applicator)
4 selectable fixed settings	at 16 Hz max. 120 mJ
Mode of operation	Intermittent use max. 6000shocks / 15min. break
Accuracy	± 20%
Shockwave hand-piece:	Ergonomic model with anodized aluminium case and fan cooling
Dimensions	230 mm in length, 50 mm diameter
Weight	850 g (with cable)
Service life	2,000,000 shocks (minimum) Applicators exchangeable without any tools (6 / 15 / 25 mm diameter)
Dimensions (complete with case)	L 580 mm / W 470 mm / H 250 mm
Total weight	12 kg (total with case)
Environmental conditions	
Operational environment	10 to 25 °C (50 to 77 °F); 700 hPa – 1060 hPa, 20% to 80% rel. humidity, not condensed
Storage / Transport	
Short-term	-10 to 55 °C (14 to 131 °F); 700 hPa – 1060 hPa, 20% to 80% rel. humidity, not condensed
Long-term	0 to 40 °C (32 to 104 °F); 700 hPa – 1060 hPa, 20% to 80% rel. humidity, not condensed

16 Safety and Performance standards

IEC 60601-1 : General requirements for the safety of electrical medical systems, including Annex 1, national differences for Australia, Canada and the United States.

Safety class according to IEC 60601-1 : class I type BF

CE₀₁₉₇

: This equipment complies with all requirements of the Medical Device Directive (93/42/EEC).

Medical device classification : IIa

17 EMC details

Medical electrical devices such as the Endopuls 811 are subject to special precautions with regard to electromagnetic compatibility (EMC) and must be installed and commissioned in accordance with the EMC advice given in the instructions for use and accompanying documents.

Portable and mobile RF communication systems (e.g. mobile phones) may interfere with medical electrical equipment.

The Endopuls 811 should only be operated with the original mains cable specified in the list of contents delivered.

Operating the device with any other mains cable can lead to increased emissions or reduced interference immunity of the device.

Guidelines and manufacturer's declaration – electromagnetic interference		
The device Endopuls 811 is intended for operation in an electromagnetic environment as indicated below. The customer or user of the Endopuls 811 should ensure that it is operated in such an environment.		
Interference tests	Conformity	Electromagnetic environment guideline
RF emissions according to CISPR 11	Group 1	The device Endopuls 811 uses RF energy solely for its internal functioning. Its RF emission is therefore very low and it is unlikely that this will cause interference to neighbouring electronic equipment.
RF emissions according to CISPR 11	Class A	The device Endopuls 811 is suitable for use in all installations including those in a residential environment and those which are directly connected to the public mains network which also supplies buildings which are used for residential purposes.
Harmonic emissions according to IEC 61000-3-2	Class A	
Voltage fluctuation emissions and flicker according to IEC 61000-3-3	Conforms	

The device should not be used when placed immediately next to or stacked on top of other devices. If operation is necessary when immediately next to or stacked on top of other devices, the device should be monitored to ensure it is operating as intended in this arrangement.

Guidance and manufacturer's declaration – Electromagnetic immunity

The Endopuls 811 device is intended for use in the electromagnetic environment specified below. The customer or the user of the Endopuls 811 device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge (ESD) to IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 6100-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T for 0.5 cycle) 40% U_T (60% dip in U_T for 5 cycles) 70% U_T (30% dip in U_T for 25 cycles) <5% U_T (>95% dip in U_T for 5 seconds)	<5% U_T (>95% dip in U_T for 0.5 cycle) 40% U_T (60% dip in U_T for 5 cycles) 70% U_T (30% dip in U_T for 25 cycles) <5% U_T (>95% dip in U_T for 5 seconds)	Mains power quality should be that of a typical commercial or hospital environment. The user of the Endopuls 811 requires continued operation during power mains interruptions. It is recommended that the Endopuls 811 be powered from an uninterruptable power supply or a battery.
Power frequency (50/60 Hz) magnetic field to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the AC mains voltage prior to application of the test level.

The main features of the Endopuls 811 are as follows: interference-free delivery of shockwaves, interference-free control of all functions. Uninterrupted operation is not required with the use intended.

Guidelines and manufacturer's declaration – electromagnetic interference immunity

The device Endopuls 811 is intended for operation in the electromagnetic environment specified below. The customer or user of the Endopuls 811 should ensure that it is used in such an environment.

Interference immunity tests	IEC 60601-test level	Compliance level	Electromagnetic environment - guidelines
<p>Conducted RF disturbance variables according to IEC 61000-4-6</p> <p>Radiated RF disturbance variables according to IEC 61000-4-3</p>	<p>3 V_{effektive value} 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V_{effektive value} 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>Portable and mobile radio equipment should not be used any closer to the Endopuls 811, including cables, than the recommended separation distance calculated from the equation applicable to the transmission frequency.</p> <p>Recommended separation distance:</p> <p>$d = 1.2 \sqrt{P}$</p> <p>$d = 0.35 \sqrt{P}$ for 80 MHz to 800 MHz</p> <p>$d = 0,7 \sqrt{P}$ for 800 MHz to 2.5 GHz</p> <p>Where P is the rated power of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>According to an investigation in situ^a, the field strength of stationary radio transmitters should be less than the compliance level at all frequencies.</p> <p>Interference may occur in the vicinity of equipment which is marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz the higher frequency range is applicable.
 NOTE 2 These guidelines may not be applicable in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Theoretically, it is not possible to exactly predict the field strengths of fixed transmitters such as base stations for radio telephones and land mobile radios, amateur radio stations, AM and FM radio and TV broadcasting. To determine the electromagnetic environment in relation to the fixed transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location where the Endopuls 811 device is to be used exceeds the above compliance levels, the Endopuls 811 device should be monitored in order to ensure that it is functioning as intended. If unusual features are noticed, additional measures may be necessary such as re-orienting or relocating the Endopuls 811 device.

b Above the frequency range from 150 kHz to 80 MHz the field strength should be less than 3 V/m.

Recommended separation distances between portable and mobile RF telecommunications equipment and the Endopuls 811 device

The Endopuls 811 device is intended for operation in an electromagnetic environment where RF disturbances are monitored. The customer or user of the Endopuls 811 device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF telecommunications equipment (transmitters) and the Endopuls 811 device – according to the output power of the communications device, as indicated below.

Rated output of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d= 1.2 \sqrt{P}$	80 MHz to 800 MHz $d= 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d= 0.7 \sqrt{P}$
0.01	0.12	0.035	0.07
0.1	0.38	0.11	0.22
1	1.2	0.35	0.70
10	3.8	1.1	2.2
100	12	3.5	7

For transmitters rated at a maximum output which is not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the respective column, whereby P is the maximum rated output of the transmitter in Watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz the higher frequency range is applicable.

NOTE 2 These guidelines may not be applicable in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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