# Sonopuls 190

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo  
Guía de usuario

Copyright:

Enraf-Nonius B.V.  
P.O. Box 12080  
Vareseweg 127  
3004 GB Rotterdam  
The Netherlands  
Tel: +31 (0)10 – 20 30 600  
Fax: +31 (0)10 – 20 30 699  
info@enraf-nonius.nl  
[www.enraf-nonius.com](http://www.enraf-nonius.com)





Part number: 1631750 - 41  
November 22, 2011

# Sonopuls 190

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo  
Guía de usuario

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## TABLE OF CONTENTS

<b>1</b>	<b>PRODUCT DESCRIPTION .....</b>	<b>33</b>
<b>2</b>	<b>FOREWORD .....</b>	<b>33</b>
2.1	This manual .....	33
2.2	Product liability .....	33
<b>3</b>	<b>PRECAUTIONARY INSTRUCTIONS .....</b>	<b>34</b>
<b>4</b>	<b>INTENDED USE ULTRASOUND THERAPY .....</b>	<b>35</b>
4.1	Indications/Contraindications and Adverse Effects for Ultrasound .....	35
4.2	Parameters .....	36
<b>5</b>	<b>PACKAGE CONTENTS .....</b>	<b>36</b>
5.1	Standard Accessories ultrasound .....	36
5.2	Optional Accessories .....	36
<b>6</b>	<b>INSTALLATION .....</b>	<b>37</b>
6.1	Connection to mains supply .....	37
<b>7</b>	<b>APPLICATION NOTES .....</b>	<b>37</b>
7.1	Ultrasound .....	37
7.1.1	Contact Control .....	37
7.1.2	The contact medium .....	37
7.1.3	Before treatment .....	37
7.1.4	During treatment .....	37
7.1.5	After treatment .....	38
7.2	Operator Controls .....	38
<b>8</b>	<b>BASIC OPERATION .....</b>	<b>40</b>
8.1	Turn on the apparatus .....	40
8.2	Home menu .....	40
8.3	System Settings .....	40
8.3.1	Set Language .....	40
8.3.2	Set Contrast .....	41
8.3.3	Select Ultrasound Norm .....	41
8.3.4	System Info .....	41

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# 1 Product Description

---

The Sonopuls 190 is a multi-frequency ultrasound unit. The multi-frequency ultrasound applicators for the 100-serie provide both 1 and 3 MHz operation. The applicators can be used with both, large and small treatment surfaces and are suitable for subaquatic treatments.

Contact control suspends the application of ultrasonic energy when acoustical contact with the treatment area becomes insufficient.

The user can connect two ultrasound applicators. These can be activated from the ultrasound menu.

**WARNING (USA only):** *The Sonopuls 190 is a prescription device that should only be used under the supervision or by the order of a physician or other licensed healthcare provider.*

---

## 2 Foreword

---

### 2.1 This manual

This manual has been written for the owners and operators of the Sonopuls 190. It contains general instructions on operation, precautionary practices, maintenance and parts information. In order to maximize the use, efficiency and lifespan of your unit, please read this manual thoroughly and become familiar with the controls as well as the accessories before operating the unit.

Specifications put forth in this manual were in effect at the time of publication. However, owing to Enraf-Nonius BV's policy of continual improvement, changes to these specifications may be made at any time without obligation on the part of Enraf-Nonius BV.

### 2.2 Product liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortcomings of the product.

To the maximum extent permitted by applicable law, in no event will Enraf-Nonius or its suppliers or resellers be liable for any indirect, special, incidental or consequential damages arising from the use of or inability to use the product, including, without limitation, damages for loss of goodwill, work and productivity, computer failure or malfunction, or any and all other commercial damages or losses, even if advised of the possibility thereof, and regardless of the legal or equitable theory (contract, tort or otherwise) upon which the claim is based. In any case, Enraf-Nonius's entire liability under any provision of this agreement shall not exceed in the aggregate the sum of the fees paid for this product and fees for support of the product received by Enraf-Nonius under a separate support agreement (if any), with the exception of death or personal injury caused by the negligence of Enraf-Nonius to the extent applicable law prohibits the limitation of damages in such cases.

The opposing party (product's user or its representative) shall disclaim Enraf-Nonius from all claims arising from third parties, whatever nature or whatever relationship to the opposing party.

### 3 *Precautionary Instructions*

---

In this section general Warnings and Precautions are listed, that you should be aware of when using the Sonopuls 190. See also chapter 4 for Warnings and Precautions that are application specific.

#### **WARNING:**

- Federal law (USA only) restricts this device to sale by, or on the order of, a physician or licensed practitioner. This device should be used only under the continued supervision of a physician or licensed practitioner.
- Make certain that the unit is electrically grounded by connecting only to a grounded electrical service receptacle conforming to the applicable national and local electrical codes.
- Do not operate the unit in an environment of short-wave or micro-wave diathermy use as this can result in burns beneath the electrodes.
- Care must be taken when operating this equipment around other equipment. Potential electromagnetic or other interference could occur to this or to the other equipment. Try to minimize this interference by not using other equipment in conjunction with it.
- This equipment is not suitable for use in the presence of flammable anesthetics mixture with air, oxygen, or nitrous oxide.
- This device should be kept out of the reach of children.

#### **CAUTION:**

- Read, understand and practice the precautionary and operating instructions. Know the limitations and hazards associated with using any electrical stimulation device. Observe the precautionary and operational decals placed on the unit.
- Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous exposure to ultrasonic energy.
- Handle ultrasound applicator with care. Inappropriate handling of the ultrasound applicator may adversely affect its characteristics.
- Inspect ultrasound applicator for cracks which may allow the ingress of conductive fluid before each use.
- Inspect ultrasound applicator cables and associated connectors before each use.
- Do not operate the Sonopuls 190 when connected to any unit other than Enraf-Nonius BV devices.
- This unit should be operated, transported and stored in temperatures between 10 °C and 40 °C (50 °F and 104 °F), with a Relative Humidity ranging from 10%-100%.
- Do not expose the unit to direct sunlight, heat radiated from a heat radiator, excessive amounts of dust, moisture, vibrations and mechanical shocks.
- In the case of ingress of liquids, unplug the unit from the mains supply and have it checked by an authorized person (see the paragraph on technical maintenance)
- Before administering any treatment to a patient you should become acquainted with the operating procedures for each mode of treatment available, as well as the indications, contraindications, warnings and precautions. Consult other resources for additional information regarding the application of ultrasound therapy.

---

## 4 *Intended Use Ultrasound therapy*

---

Ultrasound is a mechanical energy consisting of high-frequency vibrations applied by means of an ultrasound applicator. These vibrations pass through the tissue of the body and are gradually absorbed and transformed into heat. The resulting temperature increase triggers biological changes to occur in the tissue for the relief of pain, relaxation of muscle spasms and reduction of joint contractures.

### 4.1 Indications/Contraindications and Adverse Effects for Ultrasound

#### Indications:

- Ultrasound is indicated for conditions that benefit from the application of deep heat: relief of pain, muscle spasms and joint contractures. The objective of therapeutic ultrasound in the treatment of selected medical conditions associated with the chronic and sub chronic conditions of bursitis/capsulitis, epicondylitis, ligament sprains, tendinitis, scar tissue healing and muscle strain, is to reduce pain.

#### Contraindications:

- The established contraindications to heat therapy itself.
- In an area of the body where a malignancy is known to be present.
- Over or near bone growth centers until bone growth is complete.
- Over the thoracic area if the patient is using a cardiac pacemaker.
- Over a healing fracture.
- Over ischemic tissues in individuals with vascular disease where the blood supply would be unable to follow the increase in metabolic demand and tissue necrosis might result.
- In the presence of metal implants of any type.
- Patients with sensory loss on the area to be treated.
- The gonads or to the developing fetus.
- The heart.
- The brains.
- The testicles.
- The eyes.
- Facial sinus as this exposes the eyes to the same hazards.
- Ultrasound should not be used on unconscious patients.

#### Precautions and Warnings:

- Precaution should be taken when using therapeutic ultrasound on patients with hemorrhagic diatheses.
- Ultrasound treatment presents a potential safety hazard in patients whose pain response has been decreased because of disease, previous surgery, ionizing radiation therapy, chemotherapy, or general or regional anaesthesia. It may cause burns. Do not use on insensitive areas or in the presence of poor circulation.
- Large thermal doses may result in regions of thermal aseptic necrosis which may not be apparent on inspection of the skin.
- See also chapter 3, Precautionary Instructions, for general Warnings and Precautions.

#### Relevant Hazards:

- Use of ultrasound in treating areas above the shoulders may pose relevant hazards. While it is recognized that certain specific conditions involving the eyes can and have been treated by qualified specialists, knowledge and experience to administer such treatments, such application carries with it recognized hazards of applying heat to the eyes.
- Treatment of the facial sinus exposes the eyes to the same hazards.
- Treatment of the thyroid, as well as lymph nodes in the neck, may expose the patient to as yet undetermined effects inasmuch as the safety of such treatments has not yet been established.

#### Potential Adverse Effects:

- Cataracts.
- Male sterility.

- Enhanced drug activity.
- Thermal stress.

## 4.2 Parameters

**Ultrasound Frequency**, expressed in MHz, is the frequency of the ultrasound waves. The ultrasound frequency determines the penetration depth, which has the largest value at 1 MHz. The ultrasound frequency can be set at 1 MHz or 3 MHz.

**Duty Cycle**, expressed in %, defines the ratio of the pulse duration to the pulse repetition time. Ultrasound can be applied in pulsed or in continuous mode.

**Effective Radiation Area** (ERA) expressed in cm<sup>2</sup>, defines the cross-sectional area of the ultrasound beam (See technical specifications for details). The Effective Radiation Area is fixed and defined by the size of the ultrasound applicator.

**Ultrasound Power** is the ultrasound output expressed in W. The ultrasound output display can be toggled between W and W/cm<sup>2</sup>. In pulsed mode the power during the pulse is displayed. The time averaged power can be obtained by multiplying this value with the Duty Cycle.

**Ultrasound Amplitude**, expressed in W/cm<sup>2</sup>, is the quotient of Ultrasound Power and Effective Radiation Area. The ultrasound output display can be toggled between W and W/cm<sup>2</sup>. In pulsed mode the Amplitude during the pulse is displayed. The time averaged Amplitude can be obtained by multiplying this value with the Duty Cycle.

---

## 5 Package Contents

---

- |          |   |
|----------|---|
| 1631.901 | Sonopuls 190, including a large ultrasound applicator           |
| 1631.902 | Sonopuls 190, including a small ultrasound applicator           |
| 1631.903 | Sonopuls 190, including a small and large ultrasound applicator |

### 5.1 Standard Accessories ultrasound

- |          |                                     |
|----------|-------------------------------------|
| 3444.929 | Contact-gel, bottle of 250 ml, 1/12 |
| 3444.357 | Mains connecting cable              |

### 5.2 Optional Accessories

- |         |   |
|---------|---|
| 3442929 | Flacon Echoston contact-gel, bottle 250 ml, box of 12 |
| 3442930 | Flacon Echoston contact-gel, bottle 850 ml, box of 12 |
| 3442931 | Echoston contact-gel, canister of 5 l                 |
| 3442932 | Dispenser-set for 5 l canister                        |

---

## 6 Installation

---

### 6.1 Connection to mains supply

- Insert the mains cable into socket [4] and connect it to a wall socket. Connect the applicator on socket [7] or [8].

#### CAUTION:

- Do not place the unit in a location where the power cord could be tripped over or pulled out during treatment.
- Do not attempt to use the unit if it is not properly grounded. Make certain that the unit is electrically grounded by connecting it only to a grounded electrical service receptacle that is in conformity with the applicable national and local electrical codes regarding medical environments.

---

## 7 Application Notes

---

### 7.1 Ultrasound

#### 7.1.1 Contact Control

The ultrasound applicator has a contact control function that suspends treatment when the acoustical contact with the body drops below a certain level (see paragraph 10.2). The indicator light on the applicator is turned on to signal this situation. The display status symbol will show pause (  ) and you will hear a beep. During this situation the applicator emits a small amount of energy to sense restoration of acoustical contact. You may experience this when the applicator only partially contacts the body. When contact restoration is sensed the treatment is resumed at the set Amplitude.

#### 7.1.2 The contact medium

To ensure efficient transfer of energy, a contact medium is required between the ultrasound applicator and the body. Air causes virtually total reflection of the ultrasound energy. The best medium for the transfer of ultrasound energy is a gel.

The gel should be applied to the part of the body to be treated and then spread out with the ultrasound applicator.

Never apply the gel to the ultrasound applicator. The applicator will register this as acoustical contact and may emit ultrasound energy, which could damage the applicator.

If the body surface is very irregular, making it difficult to obtain good contact between the ultrasound applicator and the body, or if direct contact must be avoided (e.g. due to pain), the affected area may be treated under water (subaquial method). The water should be degassed (by previous boiling) in order to prevent air bubbles arising on the ultrasound applicator and the body.

#### 7.1.3 Before treatment

Check the patient for contraindications. See section 4.1 for details.

Test the warmth sensibility of the treatment area.

To optimize ultrasound transmission, clean the skin of the treatment area with soap or a 70% alcohol solution.

Strong hair growth has to be shaved.

#### 7.1.4 During treatment

The ultrasound applicator has to be moved constantly, also with the semi-static method. During treatment the displayed ultrasound Amplitude can vary around the set value, caused by fluctuations in acoustical coupling.

Ask the patient regularly for his/her findings. If necessary the treatment will have to be adapted. The Amplitude can be reduced or the continuous mode can be changed to pulsed mode or vice versa. When there are signs that the ultrasound transmission is bad, add more contact gel or spread it with the applicator.


**CAUTION:**

- The ultrasound applicator is a precision instrument. Great care has been taken during the development and in production to obtain the best possible beam characteristics. Rough treatment (jarring or dropping) can adversely affect these characteristics, and must therefore be avoided.

**7.1.5 After treatment**

- Clean the skin of the patient and the ultrasound applicator with a towel or tissue. To disinfect the treatment applicator we advise to use Dismozon® (art. no. 3442438) in a solution between 0,25 and 1,5%.
- Check for the effects that can be expected (for example pain, circulation and mobility).
- Ask the patient to inform the therapist of any reactions.

**7.2 Operator Controls**


**[1] LCD**  
Colour TFT display 4.3" with touch screen

**display:**

**[2] Applicator:**  
Ultrasound applicator Multi-frequency 1 and 3 MHz

**[3] Applicator holder:**  
Applicator holder. Can be mounted on the right or left side of the unit.

**[4] Mains Switch**

To switch the unit on and off


When turning the unit off the display remains visible for a few seconds (low-energy method).

**[5] Connector for mains cable**

To disconnect the device from the power line, unplug the power cord.

**[6] Type number/warning sticker**

Provides information on the apparatus, such as type and serial number, as well as connection data such as mains voltage and maximum current consumption.

**[7] Connection Ultrasound applicator A**

Ultrasound applicator connection left side

**[8] Connection Ultrasound applicator B**

Ultrasound applicator connection right side

**⚠ CAUTION:**

- Connection of accessories other than the ones specified by the manufacturer can adversely affect the safety of the patient and correct functioning of the equipment, and is therefore not permitted. For combined applications only use Enraf-Nonius type BF equipment. The very low leakage current of this type of equipment ensures absolutely safe therapy.
- The ultrasound applicator is a precision instrument. Great care has been taken during the development and in production to obtain the best possible beam characteristics. Rough treatment (jarring or dropping) can adversely affect these characteristics, and must therefore be avoided.
- The selection for activating the output is in the ultrasound menu
- Only 1 canal can be active. You cannot work simultaneously with applicator A and B.


## 8 Basic Operation

### 8.1 Turn on the apparatus

- Switch the unit On via the [4].
- The unit starts by executing the self test.
- At the end of the self test the unit enters the Home menu and is ready for use.


### 8.2 Home menu

The Sonopuls 190 is equipped with touch screen technology. Via the menu you can access to all functions. Select the desired function by pushing an icon.


### 8.3 System Settings


#### 8.3.1 Set Language

- Push the icon "System Settings".


- Select from the list Language:
- Push icon "Language"


- Select your language via the ▲▼ icon
- Push OK to confirm and to return to the previous menu

### 8.3.2 Set Contrast

Select from the Brightness:

- Push icon "Brightness"

- **Adjust the brightness via the ▲▼ icon**  
Push OK to confirm and to return to the previous menu

### 8.3.3 Select Ultrasound Norm

Here you can select the ultrasound standard. The Sonopuls offers 2 standards, the USA (FDA) and the EU standard. Default setting is the USA standard.


To select a different standard, push the Ultrasound Norm icon and select another Norm via the ▲▼ icons. Push OK icon in the screen to confirm and to return to the previous menu.

### 8.3.4 System Info

Here you can find all information about your unit. Installed software is directly displayed on the left side of the system info icon. To view all other info's, push the System info icon and you have all other information about hardware on the screen.


Push OK icon in the screen to confirm and to return to the previous menu.

## 8.4 Load Clinical Protocols


The Sonopuls has 25 treatment suggestions for the most common treatments for ultrasound.


- Push the icon "Clinical Protocols" to activate the list.

Select in the list the treatment suggestion or:


- Push icon ⓘ for more therapy information about the protocol
- Push the ▲▼ icon to scroll in the list.


- The parameter screen shows all parameters at a glance
- Push the intensity read-out to change the intensity level.

Applicator and other parameters are changeable any time


- Adjust the intensity level via the ▲▼ icon.
- Push OK icon to return in the previous menu

- Place the applicator on the treatment area (with gel).
- Contact control indicator turns off and the treatment time start (to countdown)
- If there is insufficient contact, the indicator will light up and treatment time will interrupt.

## 8.5 Manual Operation


### 8.5.1 Ultrasound Frequency

- Push the action icon "Manual Operation" to activate this function


- Select a parameter in the screen by pushing the icon.
- Push icon [1] to change the frequency.


- Select the frequency by pushing the icon.
- 1 MHz is deep
  - 3 MHz is superficial
- Push OK to confirm and to return to the previous menu


### 8.5.2 Continuous / pulsed mode

- Push the icon [3] to select continuous or pulsed mode

- Push A or B
- A is pulsed mode
  - B is continuous mode
- Push OK to confirm and to return to the previous menu


If continuous pulse mode is selected no other parameters (Frequency or Duty Cycle) are available.

### 8.5.3 Pulse Frequency


- Push the icon [5] to select pulse frequency

- Select the pulse frequency via the ▲▼ icon.
- There are 3 frequencies available 100, 48 and 16Hz
  - Push OK to confirm and to return to the previous menu

### 8.5.4 Duty Cycle


- Push the icon [6] to select Duty Cycle


Select the duty cycle via the ▲▼ icon.

- There are 6 settings available 80, 50, 33, 20, 10 and 5%
- Push OK to confirm and to return to the previous menu

### 8.5.5 Applicator selection

- Push the icon [7] to select the applicator


Select the applicator A or B by pushing the icon. Push OK to confirm and to return to the previous menu


- A is the left connection
- B is the right connection

When no applicator is connected on the selected output beeb will occur.

### 8.5.6 Units


- Push the icon [8] to select “units”


Select the desired read-out of the ultrasound intensity. There are 2 read-outs, Watt/cm<sup>2</sup> and Watt. Push OK to confirm and to return to the previous menu

### 8.5.7 Treatment time



- Push the icon [9] to select the treatment time

Select the treatment time via the ▲▼ icon. Push OK to confirm and to return to the previous menu

### 8.5.8 Intensity setting



- Push the icon [8] to select intensity setting

Select the intensity level via the ▲▼ icon.  
Push OK to confirm and to return to the previous menu

### 8.6 Storing a protocol

Prior to the execution of a modified protocol or manual setting, it can be saved for later use under Favorites.

- Push the icon [11] to store the settings


Select a program number via the ▲▼ icon.  
Push OK to confirm and to return to the previous menu

Notes:

- Program number 0 is the default start up setting.
- You can extend the name of the program with your own terminology


## 8.6.1 Load Favorites

- Push the icon Favorites

Select a program via the ▲▼ icon.

- Place the applicator on the treatment area (with gel).
- Contact control indicator turns off and the treatment time will start (countdown starts)

If there is insufficient contact, the indicator will light up and treatment time will interrupt.

---

## 9 Maintenance and Troubleshooting

---

### 9.1 User Maintenance

#### 9.1.1 Cleaning of apparatus

To clean the unit, turn it off and unplug the power supply. Clean the unit with a damp cloth. Do not use abrasive cleaners. A small amount of mild household detergent may be used, if needed. To disinfect the unit we advise to use Dismozon® pur. in a solution between 0,25 and 1,5%

#### 9.1.2 Cleaning display panel

The display panel contains an anti-reflective coating, which needs special care when cleaned. Use a soft and dry cotton cloth or micro fibre tissue to clean the panel. To remove fingerprints or grease, use a non-abrasive glass cleaning agent. Apply a small amount of the cleaning agent to a soft cotton cloth and then carefully clean the panel.

#### CAUTION:

- Do not spray the cleaning agent directly on the glass panel.
- Do not use cleaning agents that contain strong alkalis, lye, acid, detergents with fluoride or detergents with ammonia.
- Use no liquid detergents, these can damage the apparatus.

#### 9.1.3 Ultrasound applicator

To prevent corrosion, clean and dry the contact surface immediately after use. Make sure that no ultrasound gel remains on the applicator.

To disinfect the treatment applicator we advice to use Dismozon® pur. in a solution between 0,25 and 1,5%. This product can be ordered at Enraf-Nonius under the number 3442438.

We further recommend cleaning the applicator and cable daily, using lukewarm water. Check the applicator and cable regularly for damage.

## 9.2 Troubleshooting

### 9.2.1 Error code

When the apparatus is turned on, it will first execute a self test. When an error is detected, both during the self test and during normal operation, a pop-up screen will appear on the display. When the error is displayed, all outputs will be disabled. When this situation occurs, remove all cables and switch the apparatus off and on again. When the error re-appears, stop using the device and contact your supplier.

If none of the above scenarios appear to be the problem, stop using the device and contact your supplier.

## 9.3 Technical Maintenance

On request a service manual can be made available containing: spare part list, descriptions, calibration instructions and other information which will assist the user's appropriately qualified technical personnel to repair those parts of the equipment which are designated by the manufacturer as repairable.

### CAUTION:

- Electrical safety of the device relies on a proper earthed electrical connection via the power cord. It is therefore necessary to have this connection checked annually.
- To ensure continued compliance with the 21 CFR 1050.10 standard, this unit should be adjusted and safety tested once each year. Procedures laid down in the service manual should be followed. This may be carried out by your supplier, or by another agency, authorized by the manufacturer. It is also recommended that a service history record is maintained. In some countries this is even obligatory.
- Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous exposure to ultrasonic energy.

### WARNING:


- This unit operates with mains voltages. No attempt should be made to disassemble the unit. Maintenance and repair should be carried out by authorized personnel only. The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorized persons.

### 9.3.1 Restore contact control and restore factory defaults


You can optimise the contact control yourself by means of the maintenance menu. Here you can also restore the factory defaults.

### WARNING:


- When restoring the factory defaults, previously stored values can be lost.

- Push System Settings

- Select "Maintenance" via the ▲▼ icon.


- Select the appropriate item via the ▲▼ icon.

## 9.4 End of life

The Sonopuls 190 contains materials that can be recycled and/or are noxious to the environment. Specialized companies can dismount the unit and sort out these materials. When you dispose of the unit, find out about local regulations concerning waste management.

## 10 Specifications


### 10.1.1 Continuous ultrasound

Ultrasound parameters:

Ultrasound frequency 1MHz or 3MHz

### 10.1.2 Pulsed ultrasound

Ultrasound parameters:

$f_p$  Pulse frequency, 100, 48 and 16 Hz

$t$  Pulse duration, set by duty cycle  
5 – 80 %

$t$  0.5 – 8 ms

RTPA 20 – 1.25

## 10.2 Ultrasound parameters

### Generator

Peak output Amplitude:

Duty cycle 5 – 50 %	: 0 – 3 W/cm <sup>2</sup>
Duty cycle 80 %	: 0 – 2.5 W/cm <sup>2</sup>
Duty cycle 100 %	: 0 – 2 W/cm <sup>2</sup> (continuous wave)

Peak output power for 5 cm<sup>2</sup> applicator:

Duty cycle 5 – 50 %	: 0 – 15 W
Duty cycle 80 %	: 0 – 12 W
Duty cycle 100 %	: 0 – 10 W (continuous wave)

Peak output power for 0.8 cm<sup>2</sup> applicator:

Duty cycle 5 – 50 %	: 0 – 2.4 W
Duty cycle 80 %	: 0 – 2 W
Duty cycle 100 %	: 0 – 1.6 W (continuous wave)

Output meter uncertainty : ± 20 % for any output above 10 % of maximum

Pulse frequency : 100, 48 and 16 Hz ± 1 %

Duty cycle pulsed : 5 – 80 %

Pulse duration at 100 Hz : 0.5 – 8 ms ± 10 % (set by duty cycle)

Temporal Peak to Average Ratio (RTPA): 20 – 1.25 ± 10 % (set by duty cycle)

Treatment timer : 0 - 30 min ± 0.1 min, linked to contact control

Contact control level : 65%

## 5 cm<sup>2</sup> Applicator

Ultrasound frequency:

1 MHz	: 0.98 MHz ± 2 %
3 MHz	: 3.1 MHz ± 2 %

ERA (Effective Radiation Area):

IEC 60601-2-5: 2000	: 4 cm <sup>2</sup>
21 CFR 1050.10	: 5 cm <sup>2</sup>

Beam type:

1 MHz	: Collimating
3 MHz	: Collimating

BNR (Beam Non-uniformity Ratio) : 6:1 maximum

Side radiation : 10 mW/cm<sup>2</sup> maximum

## 0.8 cm<sup>2</sup> Applicator

Ultrasound frequency:

1 MHz	: 0.98 MHz ± 2 %
3 MHz	: 3.1 MHz ± 2 %

ERA (Effective Radiation Area):

IEC 60601-2-5: 2000	: 0,6 cm <sup>2</sup>
21 CFR 1050.10	: 0,8 cm <sup>2</sup>

Beam type:

1 MHz	: Collimating
3 MHz	: Diverging

BNR (Beam Non-uniformity Ratio) : 6:1 maximum

Side radiation : 10 mW/cm<sup>2</sup> maximum

## Description of ultrasound field

The spatial distribution of the radiated field is a collimated beam (diverging for the 0.8 cm<sup>2</sup> applicator at 3 MHz) of ultrasound energy, with a decreasing Amplitude at increasing distance from the applicator surface. This field distribution applies to the radiation emitted into the equivalent of an infinite medium of distilled, degassed water at 30 °C and with line voltage variations in the range of ± 10% of the rated value. The ultrasonic beam is characterized by the Effective Radiation Area (ERA) and the Beam Non-uniformity Ratio (BNR).

The Effective Radiation Area is the cross-sectional area of the ultrasound beam. Its value depends on the ultrasound standard used:

- International: IEC 60601-2-5: 2000
- USA: 21 CFR 1050.10

The Beam Non-uniformity Ratio is the ratio of the maximum ultrasound Amplitude to the average ultrasound Amplitude, measured at the Effective Radiation Area. A low BNR value is indicative for the absence of high and potentially dangerous energy concentrations.

Amplitude modulatie	Duty cycle	Pulsduur	RTPA
16Hz	5%	3.1ms	20
16Hz	10%	6.3ms	10
16Hz	20%	12.5ms	5
16Hz	33%	20.6ms	3
16Hz	50%	31.3ms	2
16Hz	80%	50ms	1.25
16Hz	100%*	62.5ms	1
48Hz	5%	1ms	20
48Hz	10%	2.1ms	10
48Hz	20%	4.2ms	5
48Hz	33%	6.9ms	3.33
48Hz	50%	10.4ms	2
48Hz	80%	16.7ms	1.25

48Hz	100%*	20.8ms	1
100Hz	5%	0.5ms	20
100Hz	10%	1ms	10
100Hz	20%	2ms	5
100Hz	33%	3.3ms	3.33
100Hz	50%	5ms	2
100Hz	80%	8ms	1.25
100Hz	100%*	10ms	1

\* = continuous mode

### 10.3 Safety and Performance standards

IEC 60601-1 : General requirements for the safety of electrical medical systems, including Annex I, national differences for Australia, Canada and the United States.

Safety class according to IEC 60601-1 : class I type BF

IEC 60601-2-5 : Particular requirements for the safety of ultrasonic therapy equipment.

: This equipment complies with all requirements of the Medical Device Directive (93/42/EEC).

Medical device classification

: IIa

21 CFR 1050.10

: This equipment complies with all requirements of 21 CFR 1050.10, Performance Standard for Ultrasonic Therapy devices.

### 10.4 EMC details

#### 10.4.1 Guidance and manufacturer's declaration – electromagnetic immunity

The Sonopuls 190 is intended for use in the electromagnetic environment specified below. The customer or the user of the Sonopuls 190 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output Lines	± 1 kV for power supply lines not applicable	Mains power quality should be that of a typical computer room.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical computer room.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles	< 5% $U_T$ for 0,5 cycle  < 5% $U_T$ for 1 cycle	Mains power quality should be that of a typical computer room. If the user of the Sonopuls 190 requires continued operation during

	70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	70% $U_T$ for 25 cycles < 5% $U_T$ For 5 sec.	power mains interruptions, it is recommended that the Sonopuls 190 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3A / m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical computer room.

NOTE  $U_T$  is the a.c.mains voltage prior to application of the test level.

#### 10.4.2 Guidance and manufacturer's declaration – electromagnetic immunity

The Sonopuls 190 is intended for use in the electromagnetic environment specified below. The customer or the user of the Sonopuls 190 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
  Conducted RF IEC 61000-4-6  Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz  3 V / m 80 MHz to 2,5 GHz	3 V  3 V / m	Portable and mobile RF communications equipment should be used no closer to any part of the Sonopuls 190, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. <b>Recommended separation distance</b> $d = [3,5/3]\sqrt{P}$ $d = [3,5/3]\sqrt{P}$ 80 MHz to 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2,5 GHz  where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

A Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Sonopuls 190 is used exceeds the applicable RF compliance level above, the Sonopuls 190 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Sonopuls 190.

B Over the frequency range 15 kHz to 8 MHz, field strengths should be less than 3 V / m.

### 10.4.3 Recommended separation distances between portable and mobile RF communications equipment and the Sonopuls 190

The Sonopuls 190 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Sonopuls 190 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Sonopuls 190 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [ W ]	Separation distance according to frequency of transmitter [ m ]		
	150 kHz to 80 MHz $d = [3,5/V_1]\sqrt{P}$	80 MHz to 800 MHz $d = [3,5/E_1]\sqrt{P}$	800 MHz to 2,5 GHz $d = [7/E_1]\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,69	3,69	7,38
100	11,67	11,67	23,33

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

## 10.5 Technical Data

Mains voltage : 100 – 240 V ± 10%  
Mains frequency : 50 / 60 Hz

Main unit:  
Dimensions : 21 x 19 x 9 cm (w x d x h)  
Weight : 760 gr.

Environmental conditions for transport and storage:  
Temperature : -20° to +70° C  
Relative humidity : 10 tot 50 % (23°C contained in original packaging)  
: 10 tot 90 % (40°C contained in original packaging)  
Atmospheric pressure : 500 to 1060 hPa

Environmental conditions for normal use:  
Temperature : 10° to 40° C  
Relative humidity : 10 to 90 % non condensing  
Atmospheric pressure : 500 to 1060 hPa

*Technical modifications reserved*

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