
# EN-TreeP

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo

Copyright:

Enraf-Nonius B.V.  
P.O. Box 12080  
3004 GB Rotterdam  
The Netherlands  
Tel: +31 (0)10 – 20 30 600  
Fax: +31 (0)10 – 20 30 699  
[info@enraf-nonius.nl](mailto:info@enraf-nonius.nl)  
[www.enraf-nonius.com](http://www.enraf-nonius.com)

Part number: 1400755-44  
4th edition, April 2009

# EN-TreeP

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo

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## 1 Preface

We very much appreciate your purchase of the EN-TreeP Pulley. The EN-TreeP Pulley offers the possibility to practice in a very flexible way medical training therapy.

The package contains the EN-TreeP Pulley and the operation manual, in which you are reading now.

In general, medical training therapy is practiced with a set of 3 products. Besides the EN-TreeP Pulley (figure 1, page 32) these are the EN-Tree Bench (fig. 2) and the EN-Tree Train (fig. 3). This equipment is intended for training and therapy under supervision of a (para)medic.

When developing and producing the EN-TreeP Pulley, much consideration was given to reliability, safety, user comfort and durability. In order to be able to use your EN-TreeP safely, it is important that you read the instructions before use.

## 2 Installation

- The EN-TreeP Pulley must be installed in such a way that its use is not limited by other objects. Keep at least 2.50 m free around the unit
- The EN-TreeP Pulley must be attached to the wall with suitable fittings. The EN-TreeP Pulley must only be mounted on a brick or concrete wall. Position the plugs firmly in the centre of a brick in case of a brick wall. The fixation must be able to resist a horizontal pulling force of 950 N
- In case of installation of 2 EN-TreeP Pulleys, we advise you to place the centre of the EN-TreeP Pulley at 2.5 m from both sides of the corner (fig. 4)

Connection of an adapter other than the prescribed charger can adversely affect the safety of the user and the functioning of the unit, and is therefore not permitted.

Do not use the EN-treeP in "wet rooms"

## 3 Directions with regard to safety

- Always use the original Enraf-Nonius accessories that are meant for this product.
- The EN-TreeP Pulley should be operated or used by one person at ones
- Limbs must be kept clear from moving parts
- Clean the cuffs after usage
- During exercising hands must be kept clear from the pneumatic cylinder
- Never re-set weights during exercising
- Avoid the cord grating along sharp parts
- Keep spectators away from the pulley when it is in use (at least 1 meter)
- The 45 degrees Pull Down accessory is only intended for support during exercising. Never apply the full body weight to the Pull Down accessory. The maximum admissible weight for the 45 degrees Pull Down accessory is 50 kilogrammes (fig. 7)
- When using the latissimusgrip, first the load must be selected before the grip is mounted. If you change the weight in the mean time you have to hold the upper weight to prevent the grip falling down
- When used by children, supervision by a knowledgeable (para)medic is obligatory. At all times keep children away from the equipment when no supervision is at hand!
- A correct biomechanical positioning of the body is necessary
- When the EN-TreeP Pulley is not used in the correct way, injury may occur. Supervision by a knowledgeable paramedic is therefore obligatory.

### 3.1 Safety aspects

Safe execution of the procedures in this manual requires technical experience and general knowledge of safety precautions as well as general knowledge of the human movement system.

Warnings, Cautions and Notes have been used throughout the manual to bring special matters to the immediate attention of the reader.

- The **Warning** concerns danger to safety of the user or technician
- The **Caution** concerns an act that might damage the equipment
- The **Note** points out a statement deserving more emphasis than the general text, but which does not deserve a **Warning** or **Caution**

The sequence of steps in a procedure is important for the safe handling of the equipment and prevention of damage. Therefore never change the sequence of a procedure.

### 3.2 Legal aspects

This manual and the information herein are the property of Enraf-Nonius B.V. (Rotterdam, The Netherlands). Enraf-Nonius B.V. disclaims any responsibility for personal injury and/or damage to equipment caused by:

- Negligence or disregard of a **Warning** or a **Caution**
- Deviation of any of the described procedures
- Execution of activities that are not described in this manual
- Ignorance of the general safety precautions for handling tools and the use of electricity
- Usage of the device in another manner as specified in this manual.

#### 3.2.1 Warnings:

- Do not change the load settings during exercising or training: load changes are only allowed in the starting position
- For safety reasons keep sufficient distance from a unit used for training; one should not be within 1 meter distance from the outer borders of a unit while it's used for training
- If pain is experienced during use: stop the training immediately and consult professional medical help

#### 3.2.2 Caution:

- EN-TreeP is powered by low voltage currency, the pneumatic resistance is fed by compressed air. Use only Enraf-Nonius supplied or prescribed spare parts. A wrong currency or a wrong pressure can result in damage of the unit
- The manufacturer of the compressor advises frequent maintenance and the use of special oil lubricants. To prevent malfunction follow the instructions in the service manual of the compressor
- Take care to prevent the compression of electrical circuits wires and air tubes between the moving parts and the unit itself
- The maximum admissible weight for the 45 degrees Pull Down accessory is 50 kilogrammes (fig. 7)

#### 3.2.3 Note:

- In case of any malfunction of your EN-TreeP contact your Enraf-Nonius supplier or Enraf-Nonius B.V. (Rotterdam, The Netherlands). Have repairs, maintenance and installations carried out by an authorized Enraf-Nonius service technician
- Repairs or maintenance performed by a non Enraf-Nonius authorized service technician and/or the use of (spare) parts that are not originally Enraf-Nonius made, results in a direct loss of any guarantee
- The maximum resistance setting for the training equipment will be dependant on the pressure of the compressor

## 4 Handling and use

### 4.1 Handling

The training equipment in the EN Strenght line have uniform operation. For the En-TreeP line the operating panel is showed as fig. 6 on page 33.

### 4.2 Track

1. Orange ON-/OFF button.
2. Display for resistance level: This display displays the set-up resistance level during set-up and training. If a programmed Chip card with more training units is used and the last Series for this unit is finished, this display will show the next sequence number or -PC-.
3. Unit of resistance level [N]. 10 [N] ~1 [kgf].
4. “-” key for decreasing the resistance level.
5. “+” key for increasing the resistance level.
6. Next key: Pressing this key will write the data of e.g. an 1RM test to the Chip card. This function only will operate when a 1RM Chip card is loaded. This key will display the sequence numbers left in a sequential order.
7. Display for time indication: This display shows the duration of the training alternated with the duration of the pause between Series. Duration of the training is only displayed when this parameter is selected instead of number of Repetitions.
8. Unit of time (seconds).
9. Display for the Number of Series and Number of Repetitions. It counts up if time is the target. The repetition counter counts down, if “Number of Repetitions” is entered during set-up, if the counter reaches “0” an audible signal sounds. If the session is continued after reaching “0” then the counter will switch over to counting upwards.
10. “#Rep.” Unit of Repetitions and “S” number of Series.
11. “-” key will decrease the number of repetitions during set-up phase.
12. “+” key will increase the number of repetitions during the set-up phase.
13. “R” key: Pulsing this key will end the current Series and reset the Repetition counter to “0”. If more Series are programmed, the program of the training unit will enter the Pause phase before the next Series. If the Reset key is pressed during the last Series for the this unit, the display for the resistance level will show the sequence number of the next training unit on the Chip card. If pressed the “R” key when all entered sessions have finished, the resistance level display will show “-PC-” indicating that one have to return with the Chip card to the personal computer.
14. Chip card LED: if the LED lights, this indicates that the cip card is in use performing a read or write operation to/from the card.

### 4.3 Use

For the sequence of training therapy we direct you to the therapy guides and reference books published on this matter. At Enraf-Nonius or its representative a publication can be ordered titled EN-Xplore, part number 1400601

## 5 Product Liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortcomings of the product.

To the maximum extent permitted by applicable law, in no event will Enraf-Nonius or its suppliers or resellers be liable for any indirect, special, incidental or consequential damages arising from the use of or inability to use the product, including, without limitation, damages for loss of goodwill, work and productivity, computer failure or malfunction, or any and all other commercial damages or losses, even if advised of the possibility thereof, and regardless of the legal or equitable theory (contract, tort or otherwise) upon which the claim is based. In any case, Enraf-Nonius's entire liability under any provision of this agreement shall not exceed in the aggregate the sum of the fees paid for this product and fees for support of the product received by Enraf-Nonius under a separate support agreement (if any), with the exception of death or personal injury caused by the negligence of Enraf-Nonius to the extent applicable law prohibits the limitation of damages in such cases.

The opposing party (product's user or its representative) shall disclaim Enraf-Nonius from all claims arising from third parties, whatever nature or whatever relationship to the opposing party.

---

## 6 Mounting of the pull down accessory

---

Please refer to page 32, figure 5 for the mounting of the Pull down accessory.

**Warning:** The maximum admissible weight for the 45 degrees Pull Down accessory is 50 kilogrammes (fig. 7).

---

## 7 Intended use

---

The EN-TreeP is a pneumatic pulley for medical training purposes.

As thus the EN-TreeP can be used for:

- Prevention: regular use should prevent cardio-respiratory and locomotor impairments
- Rehabilitation: muscular dysbalances can be adjusted / equalized by means of regular use

---

## 8 Contra-indications

---

### 8.1 Absolute contra-indications:

- Myocardial dysfunction in an instable phase
- HNP in the (sub) acute phase (with irradiation and neurological signs)
- Severe osteoporosis
- Acute fractures
- Severe mechanical instability
- Early post-surgical joint modification

### 8.2 Relative contra-indications:

- Myocardial dysfunction in a stable phase
- HNP in the sub-acute phase (with irradiation and neurological signs)
- Mild osteoporosis
- Late post-surgical joint modification
- Tumours
- Inflammations in the involved area
- Non desired articular load in the involved area (p.e. with M. Osgood-Schlatter, bursitis, etc.)
- Non desired ligamentous load in the involved area (p.e. sprain with hydrops)
- Non desired soft tissue load in the involved area (p.e. with severe atrophy)
- Cardiovascular restrictions
- Pulmonary restrictions
- Psychological restriction
- Pain of such cause that movement would not relieve

---

## 9 Maintenance

---

- The steel undercarriage can be cleaned with a damp cloth and a (non abrasive) liquid household cleaner
- To keep the leather parts flexible, you can use leather fat
- After some time, The sliding and turning parts could run less smoothly, because of pollution and dehydration. You can re-grease the Pulley with ball bearings grease (which you can obtain in a bicycle shop). If required, the guiding bars of the weights can be greased with a little oil.

- Check the cord and the cord fixations regularly for signs of wear. In case of excessive wear and tear, the cord must be replaced and use of the EN-TreeP Pulley is not recommended until the cord has been replaced.
- In order to obtain spare parts, you please contact your supplier.
- Your EN-TreeP contains materials which can be recycled and/or are noxious for the environment. At the end of the duration of life, specialized concerns can take apart the EN-TreeP and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment.

Please ensure that you are well informed of the local rules and regulations regarding to the removal of equipment and accessories.

---

## 10 Finally

---

We hope that you will use the EN-TreeP Pulley for a long time with much satisfaction and we take the liberty of drawing your attention to the many other high quality products in our programme. If you should have any questions or wish to receive our catalogue, please contact our representative.

---

## 11 Technical data

---

Safety class	: EN 957-1+A1 Class SB
medical class	: 93/42/EEC Class IIa
Dimensions	: 218 x 39 x 35 cm (HxWxD)
Weight (w. covers)	: 60 kg
Height adjustment for rope	: 0 - 200 cm.
Load increments	: 40 to 720 N, in steps of 10 N
For 1 cord end	: display value divided in 6
For 2 cord ends	: display value divided in 3

For a proper operation of the repetition counter, the movement of 1 cord should be minimum 20 cm or 10 cm with 2 cords

*Technical modifications reserved.*

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