Endomed 481 / 482 / 484

Therapy guide

Protocols

**Protocollen
Protocols
Protokolle
Protocole
Protocolos**

Copyright:

Enraf-Nonius B.V.
P.O. Box 12080
3004 GB ROTTERDAM
The Netherlands
Tel.: (0)10 - 20 30 600
Fax: (0)10 - 20 30 699

Article number 1497.760-42
December 2005

Endomed 481 / 482 / 484

Therapy guide

Protocols

**Protocollen
Protocols
Protokolle
Protocole
Protocolos**

2

CONTENTS

Nederlands	5
English	27
Deutsch	49
français	73
Español	97

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The Endomed offers the opportunity to store frequently used settings in a memory. Single programs can be stored under program numbers 0 till 9. Sequential programs can be stored under cluster number 1. The pre-programmed programs can be retrieved as from cluster number 2.

2. PAIN RELIEF	29
2.1. SUPERFICIAL PAIN (DIA.)	29
2.2. NEUROGENIC PAIN (DIA.)	29
2.3. ACUTE PHASE (MF)	29
2.4. SUBACUTE PHASE (MF)	29
2.5. CHRONIC PHASE (MF)	30
2.6. ACUTE PHASE (TENS)	30
2.7. SUBACUTE PHASE (TENS)	30
2.8. CHRONIC PHASE (TENS)	30
2.9. DEEP PAIN: SUBACUTE PHASE (TENS)	30
Points 2.10 - 2.13 only apply to the Endomed 482 and 484	
2.10. POSTHERPETIC NEURALGIA DEX. TH5-TH7 (TENS) 	31
2.11. DEEP PAIN, ACUTE PHASE (MF)	31
2.12. DEEP PAIN, SUBACUTE PHASE (MF)	31
2.13. MYALGIA (MF)	31
Points 2.14 - 2.16 only apply to the Endomed 484	
2.14. HERPES ZOSTER DEX TH5-TH7 (TENS) 	31
2.15. CRPS type 1 (TENS) 	32
2.16. PSEUDO RADICULAR SYNDROME L5-S1 (TENS) 	32
3. CIRCULATION IMPROVEMENT	33
3.1. SUPERFICIAL CIRCULATION IMPROVEMENT (DIA.)	33
3.2. CIRCULATION IMPROVEMENT (TENS)	33
3.3. CIRCULATION IMPROVEMENT (MF)	33
Point 3.4 only applies to the Endomed 482 and 484	
3.4. INTERMITTENT CLAUDICATION (TENS)	33
4. MUSCLE RELAXATION	34
4.1. SINGULAR PROGRAMME (TENS)	34
4.2. SEQUENTIAL PROGRAMME (MF)	34
5. NERVE LESIONS	35
5.1. MUSCLE STIMULATION (FARADISM)	35
5.2. MUSCLE STIMULATION (TENS)	35
5.3. MUSCLE STIMULATION: PHASIC FIBRES (INITIAL STAGE) (TENS)	35
5.4. MUSCLE STIMULATION: PHASIC FIBRES (FINAL STAGE) (MF)	35
5.5. MUSCLE STIMULATION: TONIC FIBRES (INITIAL STAGE) (MF)	36
5.6. MUSCLE STIMULATION: TONIC FIBRES (FINAL STAGE) (TENS)	36
6. MUSCLE TRAINING	37
6.1. NEUTRAL (TENS)	37
6.2. SLOW-TWITCH FIBRES (TENS)	37
6.3. FAST-TWITCH FIBRES (TENS)	37
6.4. TOTAL TRAINING (MF,TENS)	37
6.5. POST-OPERATIVE (MF)	38
6.6. ENDURANCE AEROBIC (TENS)	38
6.7. ENDURANCE (MAX. AEROBIC) (TENS)	38
6.8. MAXIMAL CONTRACTION (TENS)	38

7. MUSCLE STRENGTHENING	39
7.1. SINGULAR PROGRAMME (TENS)	39
7.2. SEQUENTIAL PROGRAMME (MF,TENS)	39
7.3. SEQUENTIAL PROGRAMME (PHASIC FIBRES) (MF,TENS)	39
7.4. SEQUENTIAL PROGRAMME (TONIC FIBRES) (MF,TENS)	40
7.5. SINGULAR PROGRAMME (MF)	40
7.6. PROGRAMME BASED ON 'RUSSIAN STIMULATION'	40
8. INFLAMMATION	41
8.1. BURSITIS (MF)	41
8.2. EPICONDYLITIS (TENS)	41
8.3. TENDINITIS (MF)	41
Point 8.4 only applies to the Endomed 482 and 484	
8.4. CAPSULITIS (MF)	41
9. URINARY INCONTINENCE	42
9.1. STRESS INCONTINENCE (TENS)	42
9.2. URGE INCONTINENCE (TENS)	42
9.3. MIXED INCONTINENCE (TENS)	42
10. TRANSVERSAL LESION OF THE SPINAL CORD AND C.V.A.	43
10.1. TROPHIC IMPROVEMENT (TENS) 	43
Point 10.2 - 10.3 only applies to the Endomed 482 and 484	
10.2. CONTRACTURE RESTRICTION (TENS) 	43
10.3. TONE NORMALIZATION HYPERTONIA (TENS) 	43
Point 10.4 - 10.10 only applies to the Endomed 484	
10.4. TROPHIC IMPROVEMENT (TENS) 	43
10.5. SPASM REDUCTION (2) (TENS) 	43
10.6. TONE NORMALIZATION HYPOTONIA (TENS) 	44
10.7. PREPARATION TO STAND (TENS) 	44
10.8. STAND (TENS) 	44
10.9. CARDIOVASCULAR (TENS) 	44
10.10. SPASM REDUCTION (1) (TENS) 	45
11. AMYOTROFIE (only applies to the Endomed 482 and 484)	46
11.1. UPPER EXTREMITY AND TRUNK (TENS)	46
11.2. LOWER EXTREMITY	46
Point 11.3 and 11.4 only applies to the Endomed 484	
11.3. ATROPHY QUADRICEPS (TENS) 	46
11.4. ATROPHY HAMSTRING (TENS) 	46
12. REVALIDATION IN SPORTS (only applies to the Endomed 484)	47
12.1. ATROPHY QUADRICEPS (TENS) 	47
12.2. ATROPHY HAMSTRING (TENS) 	47
12.3. LESION ANTERIOR CRUCIATE LIGAMENT (1) (TENS) 	47
12.4. LESION ANTERIOR CRUCIATE LIGAMENT (2) (TENS) 	47
12.5. DEFICIENCY ABDOMINAL MUSCLES (TENS) 	48
12.6. TENDINITIS ROTATOR CUFF (TENS) 	48

 For more detailed explanation and application of the marked protocols () we refer to the therapy booklet "TENS treatments with the Endomed 484: multi-channel applications.".

2.1. SUPERFICIAL PAIN (DIA.)

1. General: Programme for the treatment of pain related disorders by using the milder current forms of diodynamic current as described by dr. Bernard.
2. Indication: Superficial, sharp and/or dull pain.
3. Treatment time: 8 minutes.
4. Remarks: Diodynamic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2mA/cm^2 of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: 2 min. DF, 3 min. LP, automatic polarity reverse 3 min. LP
6. Sensation: Tickling, prickling followed by alternately light to moderate vibrations.
7. Intensity level: Clearly to strongly perceptible, preferably just below motoric threshold.

2.2. NEUROGENIC (DIA.)

1. General: Programme for the treatment of neurogenic pain by using the CP form of diodynamic current (CP = courant modulé en Courtes Périodes) and automatic polarity reverse.
2. Indication: Neuralgia, causalgia, post-herpetic neuralgia, phantom limb pain (besides pain relief, local circulation is increased).
3. Treatment time: 6 minutes.
4. Remarks: Diodynamic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2 mA/cm^2 of effective treatment area. Preferably use thick pads/sponges (min. 1cm).
5. Programme: 3 min. CP, automatic polarity reverse, 3 min. CP
6. Sensation: Alternately light and rough vibrations.
7. Intensity level: Clearly perceptible, not painful. Preferably just below motoric threshod.

2.3. ACUTE PHASE (MF)

1. General: Programme for Medium Frequency current (interferential current, bipolar mode) using varying parameters in order to prevent accommodation.
2. Indication: Acute pain in absence of muscular hypertonia or as initial treatment for patients who are extremely sensitive.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 10 kHz, AMF 100 Hz, 3 min.
6. Sensation: IF, 10 kHz, AMF 80 Hz, frequency sweep 80 Hz, sweep cycle 12/12, 15 min.
7. Intensity level: Mild, rhythmically occurring tingling sensation.

2.4. SUBACUTE PHASE (MF)

1. General: Programme for Medium Frequency current (interferential current, bipolar mode) using varying parameters in order to prevent accommodation.
2. Indication: Acute and subacute pain with symptoms of muscular hypertonia.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 6 kHz, AMF 60 Hz, 3 min.
6. Sensation: IF, 6 kHz, AMF 30 Hz, frequency sweep 50 Hz, sweep cycle 6/6, 12 min.
7. Intensity level: Clearly perceptible, varying tingling sensation.

2.5. CHRONIC PHASE (MF)

1. General: Programme for Medium Frequency current (interferential current, bipolar mode) using varying parameters in order to prevent accommodation.
2. Indication: Chronic pain disorders.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4 kHz, AMF 20 Hz, 5 min.
6. Sensation: IF, 4 kHz, AMF 1 Hz, frequency sweep 29 Hz, sweep cycle 1/1, 15 min.
7. Intensity level: Fast intermittent, tingling sensations.
7. Intensity level: Strongly perceptible up to motoric level.

2.6. ACUTE PHASE (TENS)

1. General: Application of a comfortable form of TENS current for pain relief with frequency modulations and a slow sweep in order to prevent accommodation.
2. Indication: Acute pain disorders or syndromes.
3. Treatment time: 20 minutes (can be prolonged if necessary).
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 75 µs, 80 Hz, 5 min.
6. Sensation: BF.ASYM., 75 µs, 80 Hz, frequency sweep 20 Hz, sweep cycle 6/6, 15 min.
7. Intensity level: Slightly prickling, tingling.
7. Intensity level: Just perceptible, comfortable.

2.7. SUBACUTE PHASE (TENS)

1. General: Application of a comfortable form of TENS current for pain relief with frequency modulations and a slow sweep in order to prevent accommodation.
2. Indication: For pain relief of disorders at the transition acute-chronic.
3. Treatment time: 25 minutes (can even be prolonged if necessary)
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 200 µs, 50 Hz, 5 min.
6. Sensation: BF.ASYM., 150 µs, 5 Hz, frequency sweep 90 Hz, sweep cycle 1/30, 20 min.
7. Intensity level: Rough tingling vibration.
7. Intensity level: From just perceptible up to patients tolerance level.

2.8. CHRONIC PHASE (TENS)

1. General: Programme for pain relief at a chronic stage. First part of the programme slightly aggressive followed by a more milder phase.
2. Indication: Chronic pain disorders or syndromes.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 200 µs, 5 Hz, 10 min.
6. Sensation: BF.ASYM., 100 µs, 80 Hz, burst 2 Hz, 20 min.
7. Intensity level: Rhythmic vibrations. Muscular activity.
7. Intensity level: Strongly perceptible. From motoric level up to patients tolerance level.

2.9. DEEP PAIN: SUBACUTE PHASE (TENS)

1. General: Programme for pain relief by means of burst TENS.
2. Indication: Treatment of deeply felt pain at the transition acute-chronic.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 150 µs, 80 Hz, burst 2 Hz, 30 min.
6. Sensation: Rhythmic vibrations, well tolerable. Muscular activity.
7. Intensity level: Up to motoric level. Strongly perceptible.

2.10. POSTHERPETIC NEURALGIA DEX. TH5-TH7 (TENS) (only applies to the Endomed 482 and 484)

1. General: Application of a comfortable form of TENS current for pain relief.
2. Indication: Pain in the affected segments, which is still present 30 day after eruption of the vesicles.
3. Treatment time: 30 min.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.ASYM., 75 µs, 100 Hz, 30 min.
6. Sensation: Prickling tingling.
7. Intensity level: Up to patients tolerance level.

2.11. DEEP PAIN: ACUTE PHASE (MF) (only applies to the Endomed 482 and 484)

1. General: Programme for pain relief by means of two medium frequency alternating currents, causing interference (four-polar interferential).
2. Indication: Treatment of deeply felt pain (e.g. joint pain) at an acute stage.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4p, 10 kHz, AMF 100 Hz, frequency sweep 50 Hz, sweep cycle 6/6, 20 min.
6. Sensation: Comfortable, rhythmically alternating stimulation from mild (100Hz) to strong (150 Hz).
7. Intensity level: Just perceptible (sensoric stimulation).

2.12. DEEP PAIN: SUBACUTE PHASE (MF) (only applies to the Endomed 482 and 484)

1. General: Programme for pain relief by means of two medium frequency alternating currents, causing interference (four-polar interferential).
2. Indication: Treatment of deeply felt pain at the transition acute-chronic.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4p, 6 kHz, AMF 80 Hz, frequency sweep 40 Hz, sweep cycle 6/6, 20 min.
6. Sensation: Comfortable, rhythmically alternating stimulation.
7. Intensity level: Just perceptible (sensoric stimulation).

2.13. MYALGIA (MF) (only applies to the Endomed 482 and 484)

1. General: Programme for muscular pain relief by means of two medium frequency Alternating currents, causing interference (four-polar interferential).
2. Indication: Myalgia (muscular pain).
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4p, isoplanar vector field, 6 kHz, AMF 100 Hz, frequency sweep 80 Hz, sweep cycle 6/6, 15 min.
6. Sensation: Comfortable, rhythmically alternating stiulation.
7. Intensity level: Just perceptible (sensoric stimulation).

2.14. HERPES ZOSTER DEX TH5-TH7 (TENS) (only applies to the Endomed 484)

1. General: Programme for pain relief by means of a comfortable TENS current.
2. Indication: Herpes Zoster dex TH5-TH7.
3. Treatment time: 30 min.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal. Four channels are cross wise positioned above and underneath the affected segments.
5. Programme: BF.ASYM., 20 µs, 100 Hz, 30 min.
6. Sensation: Slightly prickling, tingling.
7. Intensity level: Just perceptible, comfortable.

2.15. CRPS type 1 (TENS) (only applies to the Endomed 484)

1. General: Programme for pain relief, reduction of orthosympathetic reflexactivity and normalisation of muscle activity en mobility by means of burst TENS and concentional TENS current.
2. Indication: CRPS after Radius fracture.
3. Treatment time: 30 min.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal. The electrodes are positioned paravertebral C4-C7, sub-acromial and in the course of the n. Medianus and n. Radialis.
5. Programme: Channel 1 and 2: BF.SYM., 400 µs, 100 Hz, burst 4 Hz, 30 min.
Channel 3 and 4: BF.ASYM., 20 µs, 100 Hz, 30 min.
6. Sensation: Rhythmic vibrations. Muscular activity.
7. Intensity level: Channel 1 and 2: Tolerance level.
Channel 3 and 4: Just perceptible, comfortable.

2.16. PSEUDO RADICULAR SYNDROME L5-S1 (TENS) (only applies to the Endomed 484)

1. General: Programme for pain relief by means of TENS current.
2. Indication: Pseudo radicular syndrome L5-S1.
3. Treatment time: 20 min.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal. The electrodes of channel 1 and 2 are positioned paravertebral L5-S1, channel 3 and 4 on the dorsal side of the left and right leg.
5. Programme: Channel 1 and 2: BF.SYM., 400 µs, 50 Hz, 20 min.
Channel 3 and 4: BF.ASYM., 50 µs, 100 Hz, 20 min.
6. Sensation: Channel 1 and 2; Prickling, tingling, channel 3 and 4; Slightly prickling, tingling.
7. Intensity level: Channel 1 and 2: Tolerance level
Channel 3 and 4: Just perceptible, comfortable.

3.1. SUPERFICIAL CIRCULATION IMPROVEMENT (DIA.)

1. General: Programme for improvement of the blood flow in superficial tissues.
2. Indication: Circulatory disorders.
3. Treatment time: 6 minutes.
4. Remarks: Diadynamic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2 mA/cm² of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: 2 min. DF, 2 min. CP, automatic polarity reverse, 2 min. CP
6. Sensation: Slightly prickling, followed by alternating light to moderate vibrations.
7. Intensity level: Strongly perceptible.

3.2. CIRCULATION IMPROVEMENT (TENS)

1. General: Programme for improvement of the blood flow by means of activation of the muscle pump function.
2. Indication: Circulatory disorders.
3. Treatment time: 10 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 150 µs, 15 Hz, contraction rhythm 1/1/1/3, 10 min.
6. Sensation: Sequence of fast contractions (fast ramp-up and ramp-down) followed by brief intervals.
7. Intensity level: Clearly visible muscle contractions.

3.3. CIRCULATION IMPROVEMENT (MF)

1. General: Programme for improvement of the blood flow by means of activation of the muscle pump function. Due to the application of a medium frequency current, presumably a more profound action is obtained.
2. Indication: Circulatory disorders.
3. Treatment time: 10 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4 kHz, AMF 50 Hz, frequency sweep 50 Hz, sweep cycle 1/1, contraction rhythm 1/1/1/3, 10 min.
6. Sensation: Sequence of fast contractions (fast ramp-up and ramp-down) followed by brief intervals.
7. Intensity level: Clearly visible muscle contractions.

3.4. INTERMITTENT CLAUDICATION (TENS) (only applies to the Endomed 482 and 484)

1. General: Programme for improvement of the blood flow in the lower extremities in case of ischemia.
2. Indication: Intermittent claudication.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200 µs, 4 Hz, 30 min.
6. Sensation: Fast, rhythmically vibrating stimulation.
7. Intensity level: Clearly perceptible (sensoric level).

4.1. SINGULAR PROGRAMME (TENS)

1. General: Programme for muscle relaxation by reducing the muscular tonicity.
2. Indication: Muscular hypertonia.
3. Treatment time: 5 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200 µs, 10 Hz, frequency sweep 15 Hz, sweep cycle 6/6, 5 min.
6. Sensation: Strong vibration.
7. Intensity level: Until the required sensation is reached.

4.2. SEQUENTIAL PROGRAMME (MF)

1. General: Multiple programme for muscle relaxation.
2. Indication: Muscular hypertonia.
3. Treatment time: 6 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme:
IF, 2,5 kHz, AMF 10 Hz, 1,5 min.
IF, 2,5 kHz, AMF 10 Hz, frequency sweep 15 Hz, sweep cycle 12/12, 1,5 min.
IF, 2,5 kHz, AMF 10 Hz, frequency sweep 30 Hz, sweep cycle 12/12, 1,5 min.
IF, 2,5 kHz, AMF 10 Hz, 1,5 min.
6. Sensation: Strong vibration.
7. Intensity level: Until the required sensation is reached.

5.1. MUSCLE STIMULATION (FARADISM)

1. General: Classic protocol to evoke a tetanic contraction by means of multiple electrical pulses.
2. Indication: First stage of therapeutic stimulation of striated muscles after reinnervation.
3. Treatment time: 15 minutes.
4. Remarks: It is recommended to evaluate the state of innervation by means of S-D curves. Faradic current is a Direct Current. At higher current densities DC may cause tissue damage. Intensity should not exceed 0.2 mA/cm² of effective treatment area. Preferably use thick pads/sponges (min. 1 cm).
5. Programme: Faradic current: pulse time 1 ms, pause time 20 ms, contraction rhythm of 1/1/1/6, 15 min.
6. Sensation: Moderate to fast contractions (fast ramp-up and ramp-down) followed by relatively long intervals.
7. Intensity level: Adjust intensity level until a tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

5.2. MUSCLE STIMULATION (TENS)

1. General: Programme to evoke a comfortable interrupted tetanic contraction.
2. Indication: First stage of rehabilitation of striated muscles.
3. Treatment time: 12 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200 µs, 35 Hz, contraction rhythm 3/2/1/6, 12 min.
6. Sensation: Contraction pattern characterised by a slow ramp-up and a fast ramp-down followed by intervals.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

5.3. MUSCLE STIMULATION: PHASIC FIBRES (INITIAL STAGE) (TENS)

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of TENS current.
2. Indication: First stage of rehabilitation of striated muscles with the objective to influence the twitch speed of fast-twitch (= phasic) fibres.
3. Treatment time: 10 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 100 µs, 50 Hz, contraction rhythm 1/4/1/12, 10 min.
6. Sensation: Comfortable muscle contractions in an on/off ratio of 1:3. Fast ramp-up and ramp-down followed by long intervals.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

5.4. MUSCLE STIMULATION: PHASIC FIBRES (FINAL STAGE) (MF)

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of interferential current. This programme can be regarded as a continuation of programme 5.3.
2. Indication: Advanced stage of rehabilitation of striated muscles with the objective to influence the twitch speed of fast-twitch (= phasic) fibres.
3. Treatment time: 15 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 80 Hz, contraction rhythm 2/2/1/12, 15 min.
6. Sensation: Gradually increasing stimulation (slow ramp-up and fast ramp-down) followed by a long interval.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

5.5. MUSCLE STIMULATION: TONIC FIBRES (INITIAL STAGE) (MF)

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of interferential current.
2. Indication: First stage of rehabilitation of striated muscles with the objective to influence the twitch speed of slow-twitch (= tonic) fibres.
3. Treatment time: 15 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 20 Hz, contraction rhythm 2/5/1/16, 15 min.
6. Sensation: Gradually increasing stimulation. Relatively long hold time followed by long interval.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

5.6. MUSCLE STIMULATION: TONIC FIBRES (FINAL STAGE) (TENS)

1. General: Programme to evoke a comfortable interrupted tetanic contraction by means of interferential current. This programme can be regarded as a continuation of programme 5.5.
2. Indication: Advanced stage of rehabilitation of striated muscles with the objective to influence the twitch speed of slow-twitch (= tonic) fibres.
3. Treatment time: 15 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 300 µs, 20 Hz, contraction rhythm 2/3/1/6, 15 min.
6. Sensation: Gradually increasing stimulation (slow ramp-up and fast ramp-down) followed by an interval.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

6.1. NEUTRAL (TENS)

1. General: General programme for muscle strengthening by means of biphasic symmetrical TENS current.
2. Indication: Intensive muscle training with the objective to increase muscle strength without affecting the muscle fibre composition.
3. Treatment time: 32 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 8 sequences of each 2 min. BF.SYM., 200 µs, 45 Hz, contraction rhythm 1/4/1/9, interval time of 2 minutes after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

6.2. SLOW-TWITCH FIBRES (TENS)

1. General: Specific programme for muscle strengthening by means of biphasic symmetrical TENS current.
2. Indication: Intensive muscle training with the objective to increase muscle strength of especially the slow-twitch (=tonic) fibres.
3. Treatment time: 40 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 10 sequences of each 2 min. BF.SYM., 300 µs, 30 Hz, contraction rhythm 1/4/1/6, interval time of 2 minutes after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

6.3. FAST-TWITCH FIBRES (TENS)

1. General: Specific programme for muscle strengthening by means of biphasic symmetrical TENS current.
2. Indication: Intensive muscle training with the objective to increase muscle strength of especially the fast-twitch (=phasic) fibres.
3. Treatment time: 27 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 6 sequences of each 1½ min. BF.SYM., 100 µs, 60 Hz, contraction rhythm 1/4/1/12, interval time of 3 minutes after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

6.4. TOTAL TRAINING (MF,TENS)

1. General: Combined programme for non-specific muscle strengthening by means of alternating current. The training itself is started with and terminated by a muscle relaxation programme (comparable to warming-up and cooling-down).
2. Indication: Intensive muscle training with the objective to increase muscle strength.
3. Treatment time: 35 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 10 Hz, 1,5 min.
8 sequences of each 2 min. BF.SYM., 200 µs, 45 Hz, contraction rhythm 1/4/1/9, interval time of 2 minutes after each sequence.
IF, 2,5 kHz, AMF 10 Hz, 1,5 min.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time. Strong vibrations at the start and at the end of the programme.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

6.5. POST-OPERATIVE (MF)

1. General: General programme for muscle strengthening post-operative at which no tetanic contraction occurs.
2. Indication: First stage of muscle rehabilitation with the objective to increase muscle strength without affecting the muscle fibre composition.
3. Treatment time: 7 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programma: IF, 2 kHz, AMF 20 Hz, 7 min.
6. Sensation: Perceptible contraction.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

6.6. ENDURANCE (AEROBIC) (TENS)

1. General: Specific programme for aerobic endurance training by means of biphasic symmetrical TENS current.
2. Indication: Muscle rehabilitation with the objective to increase the endurance of the muscle.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programma: BF.SYM., 200 µs, 35 Hz, contraction rhythm 1/5/0/6, 20 min.
6. Sensation: Gradually increasing stimulation followed by an interval.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

6.7. ENDURANCE (MAX. AEROBIC) (TENS)

1. General: Specific programme for maximal aerobic endurance training by means of biphasic symmetrical TENS current.
2. Indication: Muscle rehabilitation with the objective to increase the maximal endurance of the muscle.
3. Treatment time: 10 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programma: BF.SYM., 200 µs, 50 Hz, contraction rhythm 1/5/0/15, 10 min.
6. Sensation: Gradually increasing stimulation followed by an interval.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

6.8. MAXIMAL CONTRACTION (TENS)

1. General: Specific programme for training of the maximal contraction of a muscle (-group).
2. Indication: Intensive muscle training with the objective to increase the maximal contraction of a muscle.
3. Treatment time: 7 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programma: BF.SYM., 200 µs, 70 Hz, contraction rhythm 1/5/0/25, 7 min.
6. Sensation: Gradually increasing stimulation followed by an interval.
7. Intensity level: Adjust intensity level until the required contraction is reached. Do not increase intensity during operation. Discontinue treatment if contraction strength weakens.

7.1. SINGULAR PROGRAMME (TENS)

1. General: Simple programme for non-specific muscle strengthening.
2. Indication: In conditions where an increment of muscle force is required for non-injured and/or fully recovered muscles.
3. Treatment time: 33 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 6 sequences of each 150 sec. BF.SYM., 400 µs, 35 Hz, contraction rhythm 2/3/1/12, interval time of 3 minutes after each sequence.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (slow ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level until a smooth tetanic contraction occurs. Do not increase intensity during operation.

7.2. SEQUENTIAL MUSCLE STRENGTHENING PROGRAMME (MF,TENS)

1. General: Complete programme for non-specific muscle strengthening, composed of a 'warming-up', 2 training phases, multiple intervals and 'cooling-down'.
2. Indication: In conditions, where an increment of muscle force is required for non-injured and/or fully recovered muscles.
3. Treatment time: 22 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 10 Hz, followed by an interval of 1 min.
BF.SYM., 400 µs, 20 Hz, contraction rhythm 2/5/1/8, 4 min, followed by an interval of 2 min.
4 sequences of each 2 min. BF.SYM., 400 µs, 20 Hz, frequency sweep 45 Hz, sweep cycle 12/12, followed by an interval of 1 min.
IF, 2,5 kHz, AMF 10 Hz 1,5 min.
6. Sensation: Sequence of gradually increasing and decreasing stimulation, interrupted by intervals of different duration. Strong vibrations during first and final 90 seconds.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Then adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Discontinue treatment as soon as contraction strength drops.

7.3. SEQUENTIAL PROGRAMME (PHASIC FIBRES) (MF,TENS)

1. General: Complete programme for specific muscle strengthening, composed of a 'warming-up', 2 training phases, multiple intervals and 'cooling-down'.
2. Indication: In conditions, where an increment of the strength of non-injured and/or fully recovered phasic (fast-twitch) muscles is required.
3. Treatment time: 22 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2p, 2,5 kHz, AMF 10 Hz, followed by an interval of 1 min.
BF.SYM., 400 µs, 50 Hz, contraction rhythm 1/4/1/14, 4 min. followed by an interval of 2 min.
4 sequences of each 2 min. BF.SYM., 400 µs, 50 Hz, frequency sweep 100 Hz, sweep cycle 1/1, contraction rhythm 1/4/1/6, followed by an interval of 1 min.
IF, 2p, 2,5 kHz, AMF 10 Hz, 1,5 min.
6. Sensation: Sequences of relatively short, not gradually increasing/decreasing stimulation, interrupted by intervals of different duration. Strong vibrations during first and final part.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Then adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Discontinue treatment as soon as contraction strength drops.

7.4. SEQUENTIAL PROGRAMME (TONIC FIBRES) (MF,TENS)

1. General: Complete programme for specific muscle strengthening, composed of a 'warming-up', 2 training phases, multiple intervals and 'cooling-down'.
2. Indication: In conditions, where an increment of the strength of non-injured and/or fully recovered tonic (slow-twitch) muscles is required.
3. Treatment time: 22 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2p, 2,5 kHz, AMF 10 Hz, 1,5 min, followed by an interval of 1 min.
BF.SYM., 400 µs, 20 Hz, contraction rhythm of 1/6/1/8, 4 min., followed by an interval of 2 min.
4 sequences of each 2 min. BF.SYM., 400 µs, 10 Hz, frequency sweep 25 Hz, sweep cycle 12/12, followed by interval of 1 min.
IF, 2,5 kHz, AMF10 Hz, 1,5 min.
6. Sensation: Sequences of relatively short, not gradually increasing/decreasing stimulation, followed by sequences with alternating rising and falling stimulation. Interrupted by intervals of different duration. Strong vibrations during first and final part.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Than adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Discontinue treatment as soon as contraction strength drops.

7.5. SINGULAR PROGRAMME (MF)

1. General: Simple programme for non-specific muscle strengthening.
2. Indication: In conditions where an increment of muscle force is required for non-injured and/or fully recovered muscles.
3. Treatment time: 15 minutes
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: 10 sequences of each 1 min. IF, 2,5 kHz, AMF 20 Hz, frequency sweep 65 Hz, sweep cycle 1/1, contraction rhythm 1/4/1/14, followed by an interval of 30 seconds.
6. Sensation: Sequence of gradually increasing and decreasing stimulation (fast ramp-up and fast ramp-down) followed by a long interval time.
7. Intensity level: Adjust intensity level (during contraction phase) until a smooth tetanic contraction occurs. Do not increase intensity any further during operation.

7.6. PROGRAMME BASED ON 'RUSSIAN STIMULATION'

1. General: Intensive programme for muscle strengthening by means of interrupted medium frequency current. The training itself is started with and terminated by a muscle relaxation programme ('warming-up' and 'cooling-down').
2. Indication: In conditions where an increment of muscle force is required for healthy or fully recovered muscles.
3. Treatment time: 23 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2,5 kHz, AMF 10 Hz, 1,5 min.
8 sequences of each 1½ min. RS, 2,5 kHz, 50 Hz, ratio 1:1, contraction rhythm 1/7/1/9, followed by an interval of 1 minute.
IF, 2,5 kHz, AMF 10 Hz, 1,5 min.
6. Sensation: Sequences of fast, powerful contractions, interrupted by intervals. Strong vibrations at the start and at the end of the programme.
7. Intensity level: First 90 seconds: adjust intensity until a strong vibrating sensation is felt. Than adjust intensity level (during contraction phase) until patient's tolerance level. Discontinue treatment if contraction strength weakens.

8.1. BURSITIS (MF)

1. General: Programme for the improvement of the blood flow by means of medium frequency current (bipolar interferential), as part of a treatment protocol to reduce the inflammatory state.
2. Indication: Subacute and chronic bursitis, e.g. bursitis praे-patellaris.
3. Treatment time: 14 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2p, 2,5 kHz, AMF 80 Hz, frequency sweep 50 Hz, sweep cycle 6/6, 14 min.
6. Sensation: Strong, continuous stimulation.
7. Intensity level: Clearly perceptible up to strongly perceptible.

8.2. EPICONDYLITIS (TENS)

1. General: Programme for the improvement of the blood flow by means of TENS current, as part of a treatment protocol to reduce the symptoms of a 'tennis elbow'.
2. Indication: Subacute and chronic tendinitis and/or periostitis of the epicondylus, e.g. epicondylitis lateralis humeri.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 75 µs, 80 Hz, frequency sweep 20 Hz, sweep cycle 6/6, 15 min.
6. Sensation: Comfortable, fast rhythmic stimulation.
7. Intensity level: Strongly perceptible up to patient's tolerance level.

8.3. TENDINITIS (MF)

1. General: Programme for the improvement of the blood flow by means of a medium frequency current (bipolar interferential), as part of a treatment protocol to reduce the inflammatory state of a tendon.
2. Indication: Subacute and chronic tendinitis.
3. Treatment time: 14 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 2p, 2,5 Hz, AMF 80 Hz, frequency sweep 80 Hz, sweep cycle 1/30, 14 min.
6. Sensation: Intruding, strong vibrations.
7. Intensity level: Clearly perceptible up to strongly perceptible.

8.4. CAPSULITIS (MF) (only applies to the Endomed 482 and 484)

1. General: Programme for the improvement of the blood flow by means of medium frequency currents (four polar interferential), as part of a treatment protocol to reduce the inflammatory state of the articular capsule.
2. Indication: Subacute and chronic capsulitis, e.g. capsulitis of the shoulder joint.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: IF, 4p, 4 kHz, AMF 100 Hz, frequency sweep 80 Hz, sweep cycle 6/6, 15 min.
6. Sensation: Tolerable, fast rhythmic stimulation.
7. Intensity level: Clearly perceptible up to strongly perceptible.

9.1. STRESS INCONTINENCE (TENS)

1. General: Programme for improving the strength of the striated pelvic floor muscles by means of a comfortable TENS current.
2. Indication: Stress incontinence.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 150 µs, 35 Hz, contraction rhythm 0/7/0/25, 15 min.
6. Sensation: Comfortable, but clearly perceptible contraction of the pelvic floor, followed by relaxation.
7. Intensity level: Adjust intensity level (during contraction phase) until a smooth, clearly perceptible contraction of the pelvic floor occurs. The intensity is not further increased during operation.

9.2. URGE INCONTINENCE (TENS)

1. General: Programme for the inhibition of the bladder contraction and/or normalization of the tone by means of a comfortable TENS current.
2. Indication: Urge incontinence.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 100 µs, 5 Hz, 15 min.
6. Sensation: Comfortable tingling in a low beat frequency.
7. Intensity level: Adjust the intensity level until clearly perceptible (sensoric stimulation).

9.3. MIXED INCONTINENCE (TENS)

1. General: Programme, in which successively pelvic floor muscles are activated and the bladder contraction is inhibited.
2. Indication: Mixed incontinence.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 150 µs, 35 Hz, contraction rhythm 0/7/0/25, 10 min.
BF.SYM., 100 µs, 5 Hz, contraction rhythm 0/7/0/25, 10 min.
6. Sensation: Clearly perceptible contraction followed by a comfortable tingling sensation.
7. Intensity level: First adjust the intensity level a clearly perceptible contraction of the pelvic floor occurs (motoric level). The intensity during the second phase should be set at clearly perceptible (sensoric stimulation).

10.1. TROPHIC IMPROVEMENT (TENS)

1. General: Programme to improve the trophic of the tissues in the shoulder.
2. Indication: Subluxation (and pain) of the shoulder due to a C.V.A.
3. Treatment time: 90 minutes (the training must be built up from 1,5 to 6 hours a day).
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 250 µs, 35 Hz, contraction rhythm 1/10/0/12, 90 min.
6. Sensation: Increasing stimulation followed by a long interval.
7. Intensity level: Adjust intensity level until an elevation in combination with an abduction and extension of the Humerus occurs.

10.2. CONTRACTURE RESTRICTION (TENS) (only applies to the Endomed 482 and 484)

1. General: Programme for the restriction of contractures by means of stimulation of the antagonist.
2. Indication: Flexion contracture of the forearm due to a C.V.A.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 250 µs, 35 Hz, contraction rhythm 1/4/1/0, 30 min.
6. Sensation: Short increasing and decreasing stimulation.
7. Intensity level: Adjust intensity level until a clearly perceptible contraction occurs which must be pain free.

10.3. TONE NORMALIZATION HYPERTONIA (TENS) (only applies to the Endomed 482 and 484)

1. General: Programme for tone normalization by means of stimulation of the extensors of the forearm and hand for function improvement.
2. Indication: Hypertonia in the flexors of the forearm due to a C.V.A.
3. Treatment time: 10 minutes (the training must be built up from 10 to 20 min.).
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 250 µs, 15 Hz, contraction rhythm 1/8/1/0, 10 min.
6. Sensation: Increasing and decreasing stimulation.
7. Intensity level: Adjust intensity level until a clearly perceptible contraction occurs which must be pain free.

10.4. TROPHIC IMPROVEMENT (TENS) (only applies to the Endomed 484)

1. General: Programme to improve the trophic by means of circulation improvement in the muscles and surrounding tissues of the forearm.
2. Indication: In-activity due to transversal lesion of the spinal cord.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal. Only channel 1 and 2 are used on the Endomed 484.
5. Programme: BF.SYM., 300 µs, 30 Hz, contraction rhythm 1/2/1/0, max. 30 min.
2 channels; alternating stimulation channel 1 and 2. For each channel alternation of a stimulation period of 4 sec. with an interval of 8 sec.
6. Sensation: Moderate to fast contractions (fast ramp-up and ramp-down) followed by relatively long intervals.
7. Intensity level: Adjust intensity level until a clear flexion and extention movement in the fingers occurs.

10.5. SPASM REDUCTION (2) (TENS) (only applies to the Endomed 484)

1. General: Programme for the reduction of spasm by stimulation of the antagonist or flexionreflex.
2. Indication: Tonic contraction in the calf muscle due to transversal lesion of the spinal cord.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal. Only channel 1 and 2 are used on the Endomed 484.
5. Programme: BF.SYM., 300 µs, 30 Hz, contraction rhythm 1/2/1/0, 30 min.
2 channels; alternating stimulation of the left and right leg. For each channel alternation of a stimulation period of 4 sec. with an interval of 8 sec.
6. Sensation: Moderate to fast contractions (fast ramp-up and ramp-down) followed by relatively long intervals.
7. Intensity level: Adjust intensity level until a flexion movement (or intention to flexion) occurs in the stimulated leg. Increase the intensity in case the flexion-reflex reduces.

10.6. TONE NORMALIZATION HYPOTONIA (TENS) (only applies to the Endomed 484)

1. General: Programme for tone normalization by means of stimulation of the extensors and flexors of the forearm and hand for function improvement.
 2. Indication: Hypotonia forearm/hand due to a C.V.A.
 3. Treatment time: 10 minutes (the training must be built up from 10 to 20 min.).
 4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal. Only channel 1, 2 and 3 are used on the Endomed 484.
 5. Programme: BF.SYM., 250 µs, 15 Hz, contraction rhythm 1/8/1/0, 10 min.
 6. Sensation: 3 channels; sequential stimulation by which first the flexors and then the extensors contract.
 7. Intensity level: Increasing and decreasing stimulation.
- Adjust intensity level until a clearly perceptible contraction occurs which must be pain free.

10.7. PREPARATION TO STAND (TENS) (only applies to the Endomed 484)

1. General: Muscle training programme of the mm. Quadriceps and Cluteus Maximus as a preparation for standing.
2. Indication: Muscle atrophy due to transversal lesion of the spinal cord.
3. Treatment time: 35 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 300 µs, 30 Hz, contraction rhythm 1/3/0/0, 5 min.
BF.SYM., 300 µs, 30 Hz, contraction rhythm 1/3/0/0, 30 min.
4 channels; alternating stimulation of the left and right leg. In each leg alternation of a stimulation period of 4 sec. with an interval of 8 sec.
6. Sensation: Moderate to fast contractions (fast ramp-up and ramp-down) followed by relatively long intervals.
7. Intensity level: Adjust the intensity level in the warming-up phase (first 5 min.) until a clear contraction, without extension of the knee, occurs. Adjust the intensity level in the training phase until the contraction is strong enough to lift the heel from the surface. If necessary increase the intensity during the treatment.

10.8. STAND (TENS) (only applies to the Endomed 484)

1. General: Programme to stand up (see 10.8).
2. Indication: Regain functionality standing for patients who suffer from transversal lesion of the spinal cord.
3. Treatment time: 60 minutes, discontinue the training in case the contraction of the mm Quadriceps reduces.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 300 µs, 30 Hz, 60 min.
4 channels; simultaneous stimulation of the left and right leg.
6. Sensation: Strong vibration.
7. Intensity level: Same intensity level as used at programme 10.8.

10.9. CARDIOVASCULAIR (TENS) (only applies to the Endomed 484)

1. General: Programme to improve the cardiovasculair endurance.
2. Indication: Inactivity due to transversal lesion of the spinal cord.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 300 µs, 30 Hz, contraction rhythm 1/1/0/0, 30 min.
4 channels; stimulation over channel 1 and 2 is alternated with stimulation over channel 3 and 4. A stimulation period of 2 sec. is followed by an interval of 4 sec.
6. Sensation: Moderate to fast contractions (fast ramp-up and ramp-down) followed by short long intervals.
7. Intensity level: Adjust the intensity level until a clear contraction is visible. The intensity may be increased during treatment to max. 140 mA. (but may not cause pain).

10.10. SPASM REDUCTION (1) (TENS) (only applies to the Endomed 484)

1. General: Programme for the reduction of spasm by cyclic stimulation of the agonist and antagonist.
2. Indication: Clonic contractions in the lower extremities due to transversal lesion of the spinal cord.
3. Treatment time: 30 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 300 µs, 30 Hz, contraction rhythm 1/2/1/0, 30 min.
4 channels; stimulation over channel 1 and 2 is alternated with stimulation over channel 3 and 4. A stimulation period of a sec. is followed by an interval of 8 sec.
6. Sensation: Moderate to fast contractions (fast ramp-up and ramp-down) followed by relatively long intervals.
7. Intensity level: Adjust the intensity level until a clear contraction in the calf muscles occurs without movement of the leg. An adequate reaction of the muscles at the start and ending of the stimulationcycle is important. When the reactions are correct the intensity may be increased.

11.1. UPPER EXTREMITY AND TRUNK (TENS) (only applies to the Endomed 482 and 484)

1. General: Trainingsprogramme after a long period of immobilisation and in-activity, composed of a 'warming-up', training phase and 'cooling down'.
2. Indication: Muscle rehabilitation with the purpose to recover the volume of the muscle(group) by means of tetanic contractions.
3. Treatment time: 45 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200 µs, 8 Hz, 10 min.
BF.SYM., 200 µs, 45 Hz, contraction rhythm 1/4/0/15, 30 min.
BF.SYM., 200 µs, 4 Hz, 5 min.
6. Sensation: Sequence of gradually increasing and decreasing stimulation.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

11.2. LOWER EXTREMITY (TENS) (only applies to the Endomed 482 and 484)

1. General: Trainingsprogramme after a long period of immobilisation and in-activity, composed of a 'warming-up', training phase and 'cooling down'.
2. Indication: Muscle rehabilitation with the purpose to recover the volume of the muscle(group) by means of tetanic contractions.
3. Treatment time: 45 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 350 µs, 8 Hz, 10 min.
BF.SYM., 350 µs, 45 Hz, contraction rhythm 1/4/0/15, 30 min.
BF.SYM., 350 µs, 4 Hz, 5 min.
6. Sensation: Sequence of gradually increasing and decreasing stimulation.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

11.3. ATROPHY QUADRICEPS (TENS) (only applies to the Endomed 484) 

1. General: Trainingsprogramme after a long period of immobilisation and in-activity of the knee.
2. Indication: Muscle rehabilitation with the purpose to recover function and volume of the m. Quadriceps by means of tetanic contractions.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
Only channel 1, 2 and 3 are used on the Endomed 484.
5. Programme: BF.SYM., 300 µs, 50 Hz, contraction rhythm 1/5/1/5, 20 min.
3 channels; simultaneous stimulation of the m. Vastus medialis, m. Vastus lateralis and m. Rectus femoris.
6. Sensation: Sequence of relatively short increasing stimulation followed by an interval.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

11.4. ATROPHY HAMSTRING (TENS) (only applies to the Endomed 484) 

1. General: Trainingsprogramme after a long period of immobilisation and in-activity of the knee.
2. Indication: Muscle rehabilitation with the purpose to recover function and volume of the m. Hamstring by means of tetanic contractions.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 300 µs, 80 Hz, contraction rhythm 1/5/1/6, 20 min.
4 channels; simultaneous stimulation of the m. Semitendinosus, m. Biceps femoris (caput longum), m. Semimembranosus, m. Biceps femoris (caput breve).
6. Sensation: Sequence of relatively short increasing stimulation followed by an interval.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

12.1. ATROPHY QUADRICEPS (only applies to the Endomed 484) 

1. General: Muscle rehabilitation programme after a muscular lesion of the m. Quadriceps.
2. Indication: Muscle training programme to recover muscle strength by means of tetanic contractions.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal. Only channel 1, 2 and 3 are used on the Endomed 484.
5. Programme: BF.SYM., 300 µs, 50 Hz, contraction rhythm 1/5/1/5, 20 min.
3 channels; simultaneous stimulation of the m. Vastus medialis, m. Vastus lateralis and m. Rectus femoris.
6. Sensation: Sequence of relatively short increasing stimulation followed by an interval.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

12.2. ATROPHY HAMSTRING (only applies to the Endomed 484) 

1. General: Muscle rehabilitation programme after a muscular lesion of the m. Hamstring.
2. Indication: Muscle training programme to recover muscle strength by means of tetanic contractions.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 300 µs, 80 Hz, contraction rhythm 1/5/1/16, 20 min.
4 channels; simultaneous stimulation of the m. Semitendinosus, m. Biceps femoris (caput longum), m. Semimembranosus, m. Biceps femoris (caput breve).
6. Sensation: Sequence of relatively short increasing stimulation followed by an interval.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

12.3. LESION ANTERIOR CRUTIATE LIGAMENT (1) (TENS) (only applies to the Endomed 484) 

1. General: Muscle training programme for the prevention of muscular disbalance of the knee.
2. Indication: Lesion of the anterior cruciate ligament without anterior sign.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: Channel 1+2: BF.SYM., 300 µs, 50 Hz, envelop 1/5/1/0
Channel 3+4: BF.SYM., 300 µs, 80 Hz, envelop 1/5/1/0
4 channels; alternating stimulation of channel 1+2 and channel 3+4. A stimulation period of 2 sec. is followed by an interval of 2 sec.
6. Sensation: Sequence of relatively short increasing stimulation followed by an interval.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

12.4. LESION ANTERIOR CRUTIATE LIGAMENT (2) (TENS) (only applies to the Endomed 484) 

1. General: Muscle training programme for the prevention of muscular disbalance of the knee.
2. Indication: Lesion of the anterior cruciate ligament with anterior sign.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: Channel 1+2: BF.SYM., 300 µs, 50 Hz, envelop 1/5/1/0
Channel 3+4: BF.SYM., 300 µs, 80 Hz, envelop 1/5/1/0
4 channels; simultaneous stimulation over four channels. A stimulation period of 9 sec. is followed by an interval of 8 sec.
Contraction of the Hamstrings (channel 3+4), which starts first, prevents the anterior sign.
6. Sensation: Sequence increasing stimulation followed by an interval.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

12.5. DEFICIENCY ABDOMINAL MUSCLES (TENS) (only applies to the Endomed 484)

1. General: Muscle training programme to improve the performance of the abdominal muscles.
2. Indication: Function recovery of the abdominal muscles in general. Specially in rehabilitation of the lumbar lordosis causing pain in the Facet Joints.
3. Treatment time: 20 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: Channel 1+2; BF.SYM., 250 µs, 30 Hz, contraction rhythm 1/6/1/0
Channel 3+4; BF.SYM., 250 µs, 30 Hz, contraction rhythm 1/5/1/0
4 channels; alternate stimulation of channel 1+2 and channel 3+4 stimulation is followed by an interval of 6 sec.
Stimulation over channel 3+4 (cranial part) starts 1 sec. later than channel 1+2 (caudal part) to achieve a backwards rotation of the pelvis. The stimulation ends at the same time.
6. Sensation: Sequence increasing stimulation followed by an interval.
7. Intensity level: Adjust the intensity level until the desired contraction occurs. Discontinue treatment if contraction strength weakens.

12.6. TENDINITIS ROTATOR CUFF (TENS) (only applies to the Endomed 484)

1. General: Programme to improve the circulation of the rotator cuff by means of caudal translation of the humerus head.
2. Indication: Tendinitis of the rotator cuff caused by a lesion of the acromioclavicular joint.
3. Treatment time: 15 minutes.
4. Remarks: No DC component, the risk for electrochemical damage (chemical burn) is minimal.
5. Programme: BF.SYM., 200 µs, 35 Hz, contraction rhythm 3/7/3/7, 15 min.
4 channels: simultaneous stimulation over channel 1+2 (m. Pectoralis major) and channel 3+4 (m. Lattissimus dorsi). A stimulation period of 13 sec. is followed by an interval of 7 sec.
6. Sensation: Sequence of relatively long increasing stimulation followed by an interval.
7. Intensity level: Adjust the intensity level until the desired contraction strength occurs which translates the Humerus head caudal.

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