
MYOFEEEDBACK

Therapieboek
Therapy manual
Therapiebuch
Manuel de thérapie
Guia de terapia
Libro di terapia

Copyright:

Enraf-Nonius B.V.
P.O. Box 12080
3004 GB ROTTERDAM
The Netherlands
Tel: +31 (0)10 – 20 30 600
Fax: +31 (0)10 – 20 30 699
info@enraf-nonius.nl
www.enraf-nonius.com

Part number: 1496.766 - 41
December 2005

Myofeedback

Therapie boek
Therapy book
Therapiebuch
Manuel de thérapie
Guia de terapia
Libro di terapia

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TABLE OF CONTENTS

Author's note.....	23
1 Foreword	24
2 Introduction	24
2.1 What is feedback?.....	24
2.2 Biofeedback versus myofeedback	25
2.3 Myofeedback versus electromyography	25
3 Muscles and muscle function	26
3.1 Physiology of the contraction	26
3.2 The EM signal	26
3.3 From contraction to action	26
3.4 How the contraction and signal are recorded	27
3.5 Advantages of using the Myomed 932.....	27
4 Parameters	28
4.1 Treatment time	28
4.2 Adjusting the sensitivity.....	28
4.3 Threshold	28
4.4 Visual and acoustic feedback	28
4.5 Zeroing	28
4.6 Filters	28
4.7 The pressure parameter	28
5 Procedure.....	29
5.1 Continuous recording of the EM and the pressure signal.....	29
5.2 Training using a contraction and resting phase	29
5.3 Sample curve and coordination training	29
5.4 Combining the EM and pressure signals	29
5.5 Combination of EM signal and muscle-training programmes	30
5.6 The function of the analysis program.....	30
5.7 Which outputs do you use on the Myomed 932?.....	30
5.8 The Myomed as a universal stimulus device	31
6 Combining myofeedback and training equipment	31
7 Research into the efficacy of myofeedback	31
8 Specifics: the areas of application	31
8.1 Orthopaedics	31
8.1.1 <i>Scoliosis</i>	31
8.1.2 <i>Arthrosis deformans</i>	31
8.1.3 <i>RSI</i>	31
8.1.4 <i>Habitual shoulder luxation</i>	32
8.2 Sports physiotherapy	32
8.2.1 <i>Quadriceps training</i>	32
8.3 neurology	33
8.3.1 <i>Facialis paresis</i>	33
8.3.2 <i>Painful shoulder in a patient with hemiplegia</i>	33
8.4 Urology.....	33
8.4.1 <i>Stress incontinence</i>	33
8.5 Psychosomatics	33
8.5.1 <i>Tension headaches</i>	33
8.6 Paediatrics	33
8.6.1 <i>Learning movements</i>	33
8.7 Lung diseases	34
8.7.1 <i>Asthma bronchiale</i>	34
Appendix: incontinence	34

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**'What has been understood stays in the mind,
whereas what has only been learned is soon forgotten.'**

'Overzicht van de physiologie van den Mensch'. Prof. Dr. J. Jongbloed, Professor at the State University of Utrecht. 1st edition 1945, 10th edition 1967.

The book from which the above quotation is taken, 'Overzicht van de physiologie van den Mensch' ('Overview of human physiology'), was recommended to me by my physiology lecturer during my physiotherapy training. 'Ladies and gentleman', he used to say, 'it is very important that you learn how to simplify things'. In order to enable you to understand just how complex the human organism is, I shall simplify the material for you to give you a clear insight into how the various systems work.

All too often in the past, therapy books have been written as if they were textbooks. I myself have also been guilty of this, giving the reader an 'overdose of theory', with the result that he or she rapidly lays the book aside as being 'too theoretical'. We sometimes find ourselves asking, 'Was the aim of the book to help the reader implement his or her skills and knowledge via the proper application of certain forms of therapy within a particular system, or was the aim the (self-) satisfaction of the author, in order to show just how knowledgeable he is about the theory of a specific field?'

The theoretical character of some therapy books *did*, however, have a certain justification. After all, there was a need for a really good textbook of physical therapy in the narrower sense.

At the beginning of the 90s through my contacts with various teams of GPs, I became familiar with the concept of a 'skill lab'. The aim of a skill lab is to train a professional in the practical application of skills and knowledge. GPs, working under the supervision of a plastic surgeon, learned how to suture conspicuous wounds without creating superfluous repair tissue or leaving a disfiguring scar. The aim of this process was not to clone pseudo-plastic surgeons, but to enable GPs to perform sutures within the limitations of their skill and knowledge. Under the motto 'understand and practise today, do it for real tomorrow', various suturing techniques were taught using pigs' trotters. We have tried to introduce this principle into the everyday clinical and first-line practice of the physiotherapist, not just in the Netherlands but in countries all over the world.

When writing this therapy book it once again became apparent that a good, clearly written textbook on myofeedback simply wasn't available. Even so, I didn't want to write a therapy book containing an enormous amount of theoretical background. This therapy book therefore approaches the therapeutic technique of myofeedback from the *practical* perspective - as a tool in the acquisition of motor skills, in support of exercise therapy.

For those of my colleagues who want to learn more about myofeedback (in both theory and practice) I recommend following a course, because learning myofeedback is an interactive experience focused primarily on the development of knowledge and understanding. Your professional body will no doubt be able to provide you with further information.

Jan Oechies, physiotherapist
Wijk bij Duurstede, August 2002

1 Foreword

An organism can only function if its feedback systems are working well. Without feedback no life is possible.

The feedback equipment used in healthcare belongs to the class of apparatus known as 'information equipment'. In general, the activities of the various organs can be observed and assessed by measuring characteristic biopotentials. Let's take the heart as an example. Due to the spatial displacements in electrical charge that take place during the cardiac cycle, potential differences arise within (and on the surface of) the body. Both the electrical field and the potential difference between the different points on the body can provide a great deal of information about the functioning of the organ concerned - in this case the heart. The heart speaks to us because it is the 'engine of life'. Without a heart function no human life is possible.

This book contains information about the biofeedback apparatus for muscles known as the Myomed 932. The Myomed not only records the electrical activity of the muscles but also their (grip) strength. Using myofeedback provides the patient with a source of motivation and stimulus aimed at encouraging his or her deeper involvement in the therapy and the motor-learning process. After all, at a very early stage, before any muscle action is even apprehensible let alone visible, the electromyography signal shows the efforts that the patient is making. Myofeedback is also an excellent tool for the physiotherapist. After all, physiotherapeutic manipulations are based on feedback. The physiotherapist can interpret the feedback *qualitatively*, while the apparatus evaluates it *quantitatively*.

Other forms of information systems are the electroencephalogram (EEG), the electromyogram (EMG) and the electro-oculogram (EOG).

2 Introduction

2.1 What is feedback?

Feedback is essential for the functioning of an organism. This bioelectrical phenomenon provides information about the extent of the organism's functioning. Feedback is the result of an interaction between the motor function, the (central) nervous system and the sensory function.


Illustration 1: the sensomotor sensory system

There are various forms of motor functions, but they can be divided into two main groups: the *spontaneous* motor functions, based on emotion and cognition, and the *external* motor functions, which are triggered by a stimulus. Motor functions arise on the basis of movements, and movements in their turn arise on the basis of contractions. Even the isolated contractions of a muscle arise on the basis of a form of feedback: the reflex.

2.2 Biofeedback versus myofeedback

Biofeedback is the collective term for all the feedback reactions within the body. In this therapy book we shall therefore refer to bioelectrical phenomena, or biopotentials. These are all the spatial displacements in electrical charge that take place within the body: body temperature, humidity, brain functions, even the potential differences between the front and rear of the eyeball (electro-oculography). Measuring the biopotentials in muscles is known as electromyography (EMG). Myofeedback is a derivative form of the observation of the EMG signal and can be used by the physiotherapist as a functional diagnostic tool as well as in evaluation and therapy.

2.3 Myofeedback versus electromyography

As with other feedback methods, the use of the EMG signal is a method for the diagnosis and treatment of patients with neuromuscular complaints. The electromyogram is a diagnostic tool. To obtain an electromyogram, mono- or bipolar needle electrodes are attached to an individual muscle or bundle of muscles in order to measure the action potentials.

Via these electrodes the fluctuations in potential of the myofibrils in the vicinity of the electrode are then recorded at rest.

In myofeedback, surface or cavity electrodes are used to record the electrical activity in a muscle or muscle group at rest and during contraction. This method is, however, dependent on the electrical impedance of the skin. Myofeedback is used both for diagnosis and as a therapy in order to provide the patient and the therapist with information about the electrical activity in the muscles. This method provides extra feedback when the natural feedback signals of the patient are (in)sufficient to achieve the contraction. In other words, it supplements the physiological feedback. This method is especially useful in the case of muscular activity that produces a contraction or movement that is invisible to the patient. An example is to be found in the muscles of the pelvic floor and anus. Myofeedback provides extra information for the motor learning process. After all, the electromyographic signal (EM signal = the electrical resultant of the neuromuscular activity) is detected both in a contraction and in the exclusively electrical activity of a muscle when there is no discernible muscle activity.


Fig. 2 The myofeedback signal as a catalyst of the motor learning process

3 Muscles and muscle function

3.1 Physiology of the contraction

Muscle is composed of muscle fibres, which, in their turn, are made up of sarcomeres. These sarcomeres contain actin and myosin filaments. It is the latter that cause the actual contraction. The end of the nerve activates the muscle by releasing the neurotransmitter acetylcholine. As a result an electrical activation of the muscle cell occurs. A contraction mechanism arises, in which bridge connections are created between the actin and the myosin. The mechanical effect of this is a contraction, in which the muscle is 'shortened'. In myofeedback we measure the electrical activation of the muscle cell.

3.2 The EM signal

The transfer of the activity from the motor nerve to the muscle occurs at several places in the muscle. These places are called motor units. The acetylcholine mentioned above causes a depolarisation. The voltage across the membrane of the motor unit is +30 microvolts at rest and -70 microvolts after activation. This difference between the resting and activation phase is measured by the myofeedback apparatus. This is the electromyographic (EM) signal. The EM signal is stochastic (i.e. extremely random and irregular). The amplitude of the signal lies between 0 and 10 microvolts (average 0 – 0.5 microvolts) and the frequency between 0 and 500 Hz (average 50 – 150 Hz). The aim is to receive an EM signal with as much information and as little 'noise' as possible. The received signal can be disturbed by the environment, other equipment functioning in the vicinity, fluorescent lighting, shifting electrode cables, excessively large electrodes and electrodes that are not properly seated or which have been put in the wrong places.

For this reason you are recommended to:

- Clean the skin and remove any coarse hair
 - Use small disposable electrodes
 - Place the electrode on the muscle belly, not on the motor point, nor in the vicinity of the muscle-tendon transfer
 - Place the electrodes 2 or 3 finger-thicknesses from each other
 - Place the reference electrode at such a distance that at least one large joint lies between the EM electrodes and the reference electrode
 - Stick all the electrode cables onto the body with tape to prevent cables from dangling/shifting
- Any remaining interference can be suppressed by using the filter built into the Myomed 932.

3.3 From contraction to action

The sensomotor chain is the foundation of a learning nervous system (9, van Cranenbergh, 1989). Movements can only be learnt on the basis of observations. Action (motor movements) arises on the basis of concentric, eccentric, supplementary and corrective muscle contractions. All the contractions perform a function in the action (motor function). This situation can be likened to an orchestra. In this analogy, the individual instrumentalists represent the individual contractions. If all the instrumentalists play their parts properly then a good performance of the music results (= the action or motor movement). If just one instrumentalist does not play well, however, then just about the whole piece of music will sound odd.

Functional integration	<i>Autonomous phase</i>
Training	
Implementation	<i>Associative phase</i>
Sensomotor sensory chain	
Reproduction	<i>Cognitive phase</i>
Biofeedback	
Biological feedback	


Fig. 3: action and reaction: training and functionality

3.4 How the contraction and signal are recorded

The skin forms the outermost boundary of our bodies. Apart from providing protection against dehydration, infection, sunlight and other harmful influences, the skin also plays an important role in the sense of touch and in temperature regulation. The very outermost layer of the skin, the epidermis, consists partially of living and partially of dead cells. The epidermis - and especially the outermost layer of dead cells - has a much higher electrical resistance than the cells in the rest of the body. Perspiration reduces this resistance, and sebum (skin fat) increases it (just like every fatty substance that is applied to the skin). The skin resistance can be reduced by degreasing the skin or by applying a conductive gel. If there are lots of dead cells (callosities) it may be necessary to scour the skin. When measuring electrical activities in the body, electrodes are attached to measure the biopotentials. This means that the potential difference between two electrodes is being measured. The high resistance of the skin is the reason why we don't actually measure the electrical potentials in the muscle but rather the electrical activity, including the transfer impedance.

To measure the bioelectrical phenomena we use electrodes. These can either be surface electrodes or electrodes that are inserted into body cavities.

3.5 Advantages of using the Myomed 932

Myofeedback provides both the therapist and the patient with information about the (electrical) activity of the musculature. As a result the patient is able to reproduce the various contractions and movements more quickly.

Because the potentials are visually displayed, the patient can see himself being 'rewarded' for the extra effort, whether this be a muscle strengthening or a muscle relaxation. The patient also sees that selected muscle functions are being performed in the requested manner. In addition, muscular activities that do not manifest as visible movements, such as those in the pelvic floor, are also detected. The Myomed includes a number of extra options that greatly reinforce the advantages mentioned above. The EM signal can be presented graphically either as a curve or as a bar chart.

The curve shows how strong the EM signal is over a specific period. This gives the patient and the therapist an insight into the course of the electrical activity. The bar chart can be used to improve the achieved values as a unit. The 'weightlifter', which can be set in the bar chart as a representation of the EM signal, has a powerful motivating effect. The curve and the bar can also be shown in combination if desired.

The acoustic signal gives the patient the opportunity to concentrate fully on the exercises without being distracted by visual information.

In the case of pelvic floor problems, the use of pressure feedback provides extra information about the functionality of the pelvic floor musculature because it directly reflects the grip strength of this musculature.

Finally, the combination of myofeedback and muscle training programmes is an example of cognitive conditioning, because the sensory information created by the electrical stimulation can be reproduced. This muscular activity can be detected and displayed by means of the EM signal.

4 Parameters

4.1 Treatment time

Although the Myomed 932 has an adjustable treatment time, a precise setting is not essential for examination and treatment. After all, what it's all about is achieving results from the therapy. A fixed time for achieving these results can never be set. It is important, however, to note what time was required to achieve the result, in order to show the progress made in the treatment. The treatment time can be set at any value up to 60 minutes.

4.2 Adjusting the sensitivity

The sensitivity to be pre-set depends on many factors. First: the location of the electrodes and the condition of the skin at those locations. If there are many dead cells in the epithelium then a higher sensitivity must be set. This also applies to a dry or greasy skin. A larger muscle group will also require a lower sensitivity. The most important thing is to set the sensitivity so that you get the best possible signal from your Myomed.

4.3 Threshold

In order to achieve a certain objective during treatment you can motivate the patient by using a threshold. In order to trigger a (stronger) contraction the patient must raise the EM signal above the pre-set threshold. In relaxation exercises, on the other hand, the aim is to keep the EM signal *below* the pre-set threshold.

4.4 Visual and acoustic feedback

As described above, the EM signal can be displayed in two ways: as a graph or as a bar chart. An acoustic signal can be added to both these forms of presentation. This acoustic signal can be used if the patient cannot see the Myomed display because of his chosen starting position or because he or she has a visual handicap.

4.5 Zeroing

Every muscle or muscle group has a certain resting activity at the beginning of an action. This resting activity can be strongly influenced by the emotions. The 'zeroing' feature can be used to eliminate this resting activity at the start of a procedure. This makes it possible to compare the values obtained for different procedures.

4.6 Filters

In the Myomed 932 there are two types of filters:

1. If, in spite of a good electrode technique, you still have a signal without optimum information, then a filter can be activated. Activating this filter (via the pre-selection menu) will suppress the interference from the environment (fluorescent lighting, etc.).
2. In the parameter screen a filter can be activated to change the graphical display. Without the filter, just 1 point is added to the graph per second. Activating the filter produces a 'quieter' graphical display, with fewer activations.

4.7 The pressure parameter

If a vaginal or anal pressure electrode is used then the results of the measurement will be expressed as pressure values. The measurement threshold lies between 10 – 360 hPa (10 – 360 mmH₂O).

5 Procedure

5.1 Continuous recording of the EM and the pressure signal

This option is used both to diagnose and to treat. It can be used both for actions (movements) and for contractions. The term 'contraction' will be used below for convenience. The patient is first instructed how to make the contraction. The therapist then practices the contraction with the patient. The therapist looks at the bar and/or curve to see whether an EM signal can be detected. If necessary the parameters are adjusted. If 'good' information is obtained then the patient is shown the result. In strength training an attempt is then made to intensify the EM signal; if relaxation is the aim of the treatment, on the other hand, then the aim is to reduce the signal. The physiotherapist can use all the professional facilitating techniques (e.g. encouragement and offering resistance), with the exception of directly touching the muscle that is being treated.

5.2 Training using a contraction and resting phase


This option is used to train the patient in the alternate tensing or relaxation of muscles or muscle groups. This is a form of interval training. In practice it seems to make sense to make the resting phase between 1.5 times and twice as long as the contraction phase. In order to encourage the transition from training phase to functional phase shown in fig. 3, it is important to experience the difference in tonus. In practice it appears that allowing a much longer resting phase and a shorter contraction (action) phase has a beneficial effect on this process.

5.3 Sample curve and coordination training

In a coordinated movement, direction, speed and strength are all adapted to suit the aim of the movement (8, Oechies, 1994). This option enables the patient to use in a coordinated manner the muscle (group) function that is obtained. This is important for the functionality of the action. A number of sample movements are shown in which the patient has to 'ration' the muscle function. A functional movement proceeds in a coordinated fashion. This means that the direction, speed and strength of the movement are adapted to suit the aim of the movement. The sensomotor sensory chain, with proprioception, observation and precisely rationed muscle actions, is necessary for this. Because observation is essential for this process, directed training can be given with the help of the sample curves, which are arranged in order of difficulty.

5.4 Combining the EM and pressure signals

Surface electrodes and vaginal/anal pressure electrodes can be used in combination. By monitoring the electromyographic activity, excessive activity of the stomach muscles or adductors of the upper limbs can be detected. This 'excessively' high activity can mask the perceptible feedback from the muscles of the perineum. In practice this option is used almost exclusively in cases of stress incontinence. Experience shows that the stress-incontinent woman cannot generate sufficient pressure force to activate the pressure electrode; there is no sensation. According to the literature, however, there often is EM activity. By strengthening the EM signal through exercises, by functional electrostimulation and by allowing the patient to see the results via the pressure electrode a cognitive conditioning is built up.

Illustrations 4 & 5: selectivity of the actions of the pelvic floor musculature compared with the stomach musculature.

5.5 Combination of EM signal and muscle-training programmes

This option is used to monitor the effects of functional electrostimulation. When treating the musculature with FES, almost all the motor units are activated during the training phase. This gives rise to an optimum contraction feedback. During the EM phase the patient is asked to reproduce this sensation. This contraction is then monitored via the EM signal. Because both stimulation and measurement are performed via the same electrodes, compensation is almost eliminated. The literature suggests that this method is very useful in people with a pelvic floor deficiency, because they find it difficult to evaluate the effect of the stimulus in relation to the contraction.

5.6 The function of the analysis program

The analysis program can be used to perform a comparison with the results of previous treatment sessions. This option offers benefits to both the physiotherapist and the patient, as both of them can then visualise any progress in the treatment. This has an encouraging effect on the patient and increases both motivation and compliance. After all, myofeedback is a part of the overall physiotherapy treatment and will always be combined with the exercises that the patient performs (e.g. at home).

5.7 Which outputs do you use on the Myomed 932?

While using the options described above you use the outputs on the left-hand side of the Myomed (see the operating manual for the device). This also applies therefore to the combination of myofeedback and functional electrotherapy. You use the attached or inserted electrodes for both the feedback and for the stimulation.

5.8 The Myomed as a universal stimulus device

The Myomed is a complete physiotherapeutic device and offers all the options for present-day electrotherapy, along with the very extensive range of feedback functions. TENS, Interference, Träbert, Diadynamic flows, galvanisation and iontophoresis, as well as muscle strengthening programmes, are available in all their contemporary forms.

6 Combining myofeedback and training equipment

Every physiotherapist knows the risks involved in using compensatory movement or musculature when using training equipment. Training the trunk musculature using a pulley as well as (certainly) quadriceps training on a leg press are just two examples of this. This can be confirmed by recording the electrical activity of both the muscle group to be trained and the muscle group that compensates or synergises the movement.

7 Research into the efficacy of myofeedback

Myofeedback has been described extensively in the literature, and many studies have shown its usefulness. In this therapy book I want to cite the results of Berghmans L.C.M et al (3, Efficacy of biofeedback, when included with pelvic floor muscle exercise treatment for genuine stress incontinence). They compared the efficacy of exercise therapy alone with exercise therapy plus biofeedback. In both groups they saw a significant improvement after twelve treatment sessions. In the group that was treated with both exercise therapy and biofeedback, however, this result was already achieved after just six treatment sessions.

We can conclude from this that biofeedback has a favourable effect on the motor learning process.

8 Specifics: the areas of application

8.1 Orthopaedics

8.1.1 Scoliosis

The development of muscle sensation is difficult but essential. Only on the basis of a developed muscle sensation can muscles be selectively trained and the posture be corrected.

8.1.2 Arthrosis deformans

The joint deformation leads to inactivity and immobility. As a result atrophy and arthrogenic instability arise. The joint is consequently more and more severely irritated. The patient must learn to re-train the musculature, but within the pain threshold.

8.1.3 RSI

In addition to repetitive movements, stress and static loads also play an important part in the development of RSI. Severe complaints of the neck, shoulders, arms and hands can result. Relaxation exercises of the neck and shoulder muscles are often found difficult. By using the sub-threshold setting the reduction in the muscular tension can be optically displayed.

8.1.4 Habitual shoulder luxation

In this disorder a rhythmic stabilisation of the muscles of the rotator cuff is important. Coordinatively this is extremely difficult.

Procedure:

contraction-resting phase. Here the patient is asked to maintain the corrected posture during the resting phase. The result is measured during the contraction phase. For this purpose a resting phase that is twice as long as the action phase is necessary.


Sensitivity:

200 microvolts and 2 channels

8.2 Sports physiotherapy

8.2.1 Quadriceps training

Stabilisation after a lesion of the anterior cruciate ligament. Movement of the knee joint is an extremely complex set of rotatory and translatory elements.

Procedure:

continuous EMG with supra-threshold. The patient must try, through contraction, to exceed the pre-set threshold.

Sensitivity:

400 microvolts with two channels, in which the electrodes on the M. Vastus Medialis are extremely important to record the essential final stretch.

In this photo, myofeedback is shown in combination with the leg press.

Procedure:

continuous EMG with 2 channels, which record the actions both of the extensors and of the flexors of the knee.

Sensitivity:

400 microvolts and 2 channels.

In this photo, myofeedback is being used in functional exercising on a treadmill.

8.3 neurology

8.3.1 Facialis paresis

Nowhere else in the body are so many different muscles to be found on such a small surface area. Facial expressions involve a complex interplay of the muscles. Training the musculature is especially difficult, because the muscles share the common feature that on one side they are attached to bone or cartilage and on the other are connected directly to the skin.

8.3.2 Painful shoulder in a patient with hemiplegia

In the foreign physiotherapeutic literature, in particular, there are many references to the use of myofeedback. In a patient with a painful shoulder the scapulo-humeral rhythm is disturbed. Because the patient has a disturbed proprioception as a result of the hemidysfunction, myofeedback can help.

8.4 Urology

8.4.1 Stress incontinence

Motor actions of the pelvic floor musculature do not lead to visible movements. The Myomed offers various possibilities for displaying the electrical activity of the musculature. In addition to the use of the various vaginal electrodes it is also possible to place surface electrodes on the musculature of the perineum.

8.5 Psychosomatics

8.5.1 Tension headaches

Tension headaches usually present as a feeling of a band wrapped around the head or as a feeling of pressure inside it. The symptoms can be caused by hypertonia and painful neck muscles. Massage is often used but is only a symptomatic treatment. Relaxation of the musculature seems to be the correct way of dealing with it (5, Bruhn, Oleson & Melgaard, 1979)

Procedure:

contraction-resting phases.

During the resting phase the relaxation exercises are performed with the patient. During the contraction phase the results are then evaluated.

Sensitivity:

60 microvolts and 2 channels.

8.6 Paediatrics


8.6.1 Learning movements

Although child physiotherapy offers many forms of facilitation for furthering the motor learning processes, myofeedback could serve as a short-term aid.

8.7 Lung diseases

8.7.1 Asthma bronchiale

The causal substrate of this complaint cannot be treated with physiotherapy. An asthmatic patient uses auxiliary respiration muscles to ventilate and, as a result, continues to breathe 'high'. The neck and shoulder musculature, in particular, are used for this purpose. Myofeedback can be a useful aid in reducing the tension in these muscles and in learning how to ventilate efficiently (6, Abdullah, 1974)(7, Kotses and Glaus, 1981).

Procedure


continuous EMG recording with 2 channels.

Sensitivity:

100 microvolts.

Appendix: incontinence

As already described, the use of myofeedback offers benefits in cases of stress incontinence. For this purpose surface electrodes can be placed on the muscles of the perineum.

View of the underside of the perineum.

Vaginal and anal cavity electrodes can also be used. There are two kinds of vaginal electrodes:

1. The EM electrode: a cavity electrode can be used to measure the EM signal and to stimulate the muscles in conjunction with a muscle-training programme. It's important that, during the stimulation, the therapist rotate the electrode for optimum results. The EM part and the FES part are placed on one side of the electrode. Because this is a dynamic electrode technique *Constant Voltage must be set on the pre-selection menu!*
2. The pressure electrode: The pressure electrode measures the 'grip' strength of the pelvic floor muscles.

It will be clear from the above text that an EM signal is easier to measure than a pressure signal in the presence of functional disturbances of the pelvic floor.

Hygiene is very important when using cavity electrodes. These electrodes must therefore be used patient-specifically. A condom can be placed over the pressure electrode.

Proper training is absolutely essential for the use of these techniques. In many countries physiotherapists are not allowed to insert these electrodes. A practical solution to this problem is for the patient to insert them herself.


Fig. 6: action and reaction, training and functionality

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