



# PASSION FOR MOTION



*Partner for Life*

Examination and treatment couches

The equipment in this catalogue is only meant to be used by competent personnel (physiotherapy, rehabilitation or associated areas). The manufacturer cannot guarantee the suitability of the products for specific therapeutic aims.

We continually strive to improve our products so that they meet market requirements and to ensure that they are certified by international testing authorities.

The products for medical applications are given the CE logo (in accordance with the Medical Devices Directives). This logo is a guarantee of good products that have been specially developed and tested for therapeutic purposes and that are manufactured according to the strictest quality standards. Ensure that you always use products for medical applications that have this logo.

All rights reserved. For availability, please contact your local distributor.

In the production of this catalogue, environmental aspects have been taken into account.

## About Enraf-Nonius

Enraf-Nonius has specialised in the development of physiotherapy and rehabilitation products and services since 1925. Be it treatment couches, advanced physiotechanical equipment, products for active and functional rehabilitation or medical fitness, Enraf-Nonius has everything the modern motion specialist desires.

An international network of well-trained product experts places our products within reach of motion specialists worldwide. Functionality, reliability and durability are synonymous with the name Enraf-Nonius – then and now. The latest production techniques, craftsmanship and our company's continuous drive to optimise its products through customer experience confirm the image of the company with a "passion for motion".

The products and services provided by Enraf-Nonius contribute to a longer, healthier and happier life.

## Our treatment couches

As regards the Enraf-Nonius treatment couches, the key words are quality, durability and innovation. Through the years our couches have constantly been developed and modified by means of the latest technologies and by listening closely to our customers' wishes and ideas. This has resulted in unrivalled, customer-specific, safe and attractive treatment couches.

Our range consists of sturdy, simple couches; couches with individually adjustable support angles, couches with a functional design, couches for your treatment rooms and cost-effective couches.

Whichever model you choose, stability, quality and reliability are always guaranteed. Treatment couches with a stylish, modern design.

It goes without saying that every model is made of the best materials. Owing to the ergonomic design of our couches, patients and the treating specialist experience the greatest possible comfort.

Enraf-Nonius couches – the real treatment couches, unique in terms of stability, reliability and durability. Designed for the professional and for years of intensive use.

2000

2025

# Content

Examination and Treatment Couches	1
Characteristics Examination and Treatment Couches	3
Manumed Optimal 3-section	5
Manumed Optimal 2-section	7
ManuXelect 3- and 2-section	9
Manumed Osteo	11
Manumed Comfort	13
Manumed Exercise	15
Manumed Traction	16
Manumed Tilt	17
Manumed Basic	18
Colours	19
Options Examination and Treatment Couches	19
Stools	20
Rolls and cushions	20
Paper roll holders and paper rolls	21
Miscellaneous	21
Couch covers	22
Fixation accessories	23
Fixation straps	24
Accessories Manumed Traction	24
Accessories Manumed Tilt	24





# ENRAF-NONIUS

## Examination and Treatment Couches

### MANUMED OPTIMAL 3-SECTION

The most versatile treatment and examination couch.



5

### MANUMED OPTIMAL 2-SECTION

The robust, simple examination and treatment couch.



7

### MANUXELECT 3- AND 2-SECTION

The space-saving, budget-friendly couch.



9

### MANUMED OSTEO

The ideal couch for osteopathic treatments.



11

### MANUMED COMFORT

The universal treatment couch with an ideal sitting position.



13





### MANUMED EXERCISE

For safe and comfortable neurological exercise therapy.



15

### MANUMED TRACTION

Traction and treatment with the very same couch!



16

### MANUMED TILT

The comfortable and safe couch for neuromotor and cardiovascular treatments.



17

### MANUMED BASIC

High-quality basic couch for examination and treatment.



18

### PORTABLE MASSAGE COUCH

The safe, easily portable couch.



18



## Characteristics Examination and Treatment Couches



### Stability

Unsurpassed stability by use of a rigid frame with swivel bearings, which prevents the couch from moving during treatment.



### Strength

Mechanical ingenuity and many years' experience in building treatment couches are manifest in the strong frame and construction, so that even the heaviest people can lie on the Manumed.



### Reliability

Workmanship and sound materials give you the reassurance of a quality product with extremely high reliability and a long life.



### Central wheel raising mechanism

The couch is quickly and easily moved with a single push on the pedal.



### Silent and fast height adjustment

The powerful motor (or the hydraulic pump) and an optimally proportioned design, provide fast and silent height adjustment.



### Minimum height 45 cm

This makes the Manumed accessible for anyone and easy to use for manipulative treatments.



### High quality finish and welded seams

Automated production and high quality requirements guarantee a perfect finish.



### Choice of sections

Choice of two to six top sections and separate arm and leg sections, providing a suitable model for every specialisation.



### Unique upholstery with wide choice of colours

Our imitation leather upholstery contains no phthalates and is flame retardant and biocompatible. With an Enraf-Nonius couch you always have a safe couch with environment-friendly upholstery. We offer a choice of rough or smooth quality imitation leather in 18 different colours. So there is always a couch that matches your interior.






### Choice of options and accessories

Choice of traction and fixation accessories, cushions and rolls for optimising the initial position and subsequent positioning of the client.



### Breathing hole to the edge

Long breathing hole till just below the edge of the head section. This means comfort for your patient and for yourself, as you can stand very close for head/neck treatments.





#### Long breathing hole

This allows the patient to lie comfortably in any position. Particularly in combination with the flexion position, it is important that the patient takes up a correct position.



#### Ergonomic sizing

Top sections that meet the ergonomic size of men and women. Allowing the client to sit correctly, while you can take up the proper working position, without unnecessary stress on your own body.



#### Ergonomic human interface

Clever positioning and movement effect of all the handles. This means that you do not need to take up any unnecessary or incorrect postures, so that the Manumed is also extremely safe for you.



#### Comfortable armrest adjustment

The armrest has a large surface area, providing optimum relaxation in the initial position. This creates a natural and comfortable treatment position for your patient.



#### Anti-slip caps

Specially designed anti-slip caps make the couch rock-solid on any surface. The couch is therefore extremely suitable for treatments such as traction and manipulation.



#### Foot space under the couch

Generous foot space (10 cm), that can also accommodate a sun lamp or patient lift beside the couch.



#### All-round operation

The all-round operation provides electric height adjustment that can be reached unseen from any position around the couch, but which is never in the way.



#### Rounded headrest

The rounded headrest (optional) makes the head/neck region easily accessible.



#### Ergo-wings

Ideal for manipulation of the spinal column with the patient in the optimum position for you, while you avoid unnecessary stress on your own body.



#### Electrical flexion position

This brings the patient smoothly into position, without any heavy support required from the therapist. It can also be used for passive manipulation of your patient.



#### Hydraulic or electric height adjustment

Select between hydraulic height adjustment (no mains lead needed) or electric height adjustment (easy, convenient and quick).

Stability

Strength

Reliability

Central wheel raising mechanism

Silent and fast height adjustment

Minimum height 45 cm

High quality finish and welded seams

Choice of sections

Choice of colours

Choice of options and accessories

Breathing hole to the edge

Long breathing hole

Ergonomic sizing

Ergonomic human interface

Comfortable armrest adjustment

Anti-slip caps

Foot space under the couch

All-round operation

Rounded headrest

Ergo-wings

Electrical flexion position

Hydraulic or electrical height adjustment





## Manumed Optimal 3-section

The versatility of the Manumed Optimal 3-section offers you, the modern and demanding physiotherapist, the perfect treatment and examination couch. Everything is possible: lying and sitting positions, Trendelenburg and flexion position, which can be adjusted during treatment with the patient simply remaining on the couch. In short, the most universal couch available, which gives both the therapist and the patient the experience of comfortable treatment. Elegant in appearance and sturdy on its retractable castors, the Manumed Optimal 3-section is the obvious choice for the therapist who makes high demands on a professional top quality treatment couch.

- 3 (to 5) adjustable sections
- Many adjustable positions, including the flexion- and (anti-) Trendelenburg positions
- Electrically adjustable pli comes standard with electrical high/low variant
- Load capacity up to 200 kg (optional to 250 kg)
- Head section folds down completely for better access of the cervical region
- Optimally positioned handles
- Long nose hole, allowing proximal or distal movement of the patient even when lying in prone position
- Optional a triangular nose hole or without nose hole
- Wheel raising mechanism for easy movement of the couch
- Height can be adjusted hydraulically (no mains lead needed) or electrically (very convenient, smooth and quick)
- Completely relaxed starting position in the face-down position with separate head section (optional)





## Manumed Optimal 2-section

With the Manumed Optimal 2-section treatment couch you have the modern Manumed Optimal quality couch at your disposal. Robust in appearance and simple to operate. Stable and durable.

Less is more with the Manumed Optimal 2-section couch. The 2-part head section or the 2-part head section with separate arm supports offer a full range of options to the therapist looking for an uncomplicated couch for safe treatment.

- 2 (up to 4) adjustable sections
- Load capacity up to 200 kg (optional 250 kg)
- Head section folds down completely for better access of the cervical region
- Optimally positioned handles
- Long nose hole, allowing proximal or distal movement of the patient even when lying in prone position
- Optional a triangular nose hole or without nose hole
- Wheel raising mechanism for easy movement of the couch
- Height can be adjusted hydraulically (no mains lead needed) or electrically (very convenient, smooth and quick)
- Completely relaxed starting position in the face-down position with separate head section with adjustable arm rests (optional)







# MANUXELECT

## ManuXelect 3- and 2-section

The examination and treatment couch par excellence for a small treatment room and a limited budget. Because of the scissors mechanism the ManuXelect is raised fully vertically. Consequently the ManuXelect takes up only 1.4 m<sup>2</sup>! The easily adjustable sections make this couch easy to operate and stable throughout. A cost-effective variant for someone who has little space available. The ManuXelect treatment couch comes in two 2-section models and in a 3-section variant.

- Electrical height adjustment from 48.5 to 105.5 cm
- Section adjustment with gas springs
- Ergonomic handles for operation with one hand
- Scissors movement saves space. The couch travels up and down linearly
- Padded sections with soft edges
- The nose hole (with optional plug) ensures patient comfort even when lying face down









*"The complete couch, saving space and very attractively priced."*


#### MANUSELECT 3-SECTION



#### MANUSELECT 2-SECTION HEAD-SECTION



#### MANUSELECT 2-SECTION BACK-SECTION



#### TECHNICAL SPECIFICATIONS

<b>Size tabletop (lxw):</b>	207 - 209 x 67 cm
<b>Height adjustment:</b>	48,5-105,5 cm
<b>Lifting time (min.-max.):</b>	approximately 35 sec.
<b>Lifting capacity:</b>	150 kg (optional 250 kg)
<b>Power supply:</b>	100 / 120 / 230V AC 50/60 Hz
<b>Current consumption:</b>	1.5 A-350W max. (light) or 2.0 A-470W max. (heavy)
<b>Colour code frame:</b>	grey = RAL 7035



#### UPHOLSTERY COLOURS

013	004	199	169	812	852
-----	-----	-----	-----	-----	-----

(1) = Wheel raising set and plug for nosehole optional





# MANUMED

## Manumed Osteo


The ultimate in osteopathy treatment couches. With anatomically divided and individually adjustable sections, Manumed Osteo is the ideal couch for osteopathy treatments. The Manumed Osteo couch is characterised by a combination of unsurpassed stability and adjustability. The styling and the operating mechanism ensure the greatest comfort during treatment, both for the specialist and the patient. For even more operating convenience, there is also the Osteo DeLuxe, a treatment couch with no less than 3 motors with which all sections can be individually adjusted electrically.

- 4 (up to 7) independently adjustable sections, allowing an infinite number of sitting and lying positions
- Flexion and (anti-)Trendelenburg positions
- Unrivalled stability combined with adjustability
- Electrical height adjustment (and the option of sectional adjustment as well) for maximum comfort
- Completely relaxed starting position in the face-down position with separate head section (optional)
- Independently adjustable leg sections (optional)
- Padded sections with soft edges





#### MANUMED OSTEO - OPTIONAL WITH 3 MOTORS



Type	
416	Manumed Osteo, hydraulic H/L
426	Manumed Osteo, electric H/L
446	Manumed Osteo with electric section adjustment, electric H/L

#### MANUMED OSTEO, HEAD-SECTION WITH ADJUSTABLE ARM SECTIONS



Type	
418	Manumed Osteo with adjustable arm sections, hydraulic H/L
428	Manumed Osteo with adjustable arm sections, electric H/L
448	Manumed Osteo with adjustable arm sections and with electric section adjustment, electric H/L

#### MANUMED OSTEO, HEAD-SECTION WITH ADJUSTABLE ARM SECTIONS AND ADJUSTABLE LEG SECTIONS




Type	
419	Manumed Osteo with adjustable arm sections and adjustable leg sections, hydraulic H/L
429	Manumed Osteo with adjustable arm sections and adjustable leg sections, electric H/L
449	with adjustable arm sections and adjustable leg sections and with electric section adjustment, electric H/L

#### TECHNICAL SPECIFICATIONS

<b>Size tabletop (lxw):</b>	208 x 67 cm. Types 416 and 426 optional 57 or 80 cm width. Types 418 and 428 optional 80 cm width.
<b>Height adjustment:</b>	48-98 cm
<b>Lifting time (min.-max.):</b>	approximately 25 sec.
<b>Lifting capacity:</b>	150 kg
<b>Force hydraulic pump:</b>	10.000 N
<b>Force electrical motor:</b>	6.000 N
<b>Power supply:</b>	100 / 120 / 230V AC 50/60 Hz
<b>Current consumption:</b>	1.0 A max..




#### FRAME COLOURS



#### UPHOLSTERY COLOURS

012	013	004	169	141	199
155	138	019	139	127	135
811	812	852	888	933	851



# MANUMED

## Manumed Comfort

The treatment couch for everyone. The large back section makes this couch a combination of adjustability and sitting comfort. You can change from lying face up to sitting without any need for the patient to get off. This makes this couch popular not only with physiotherapists, but also with medical practitioners and aestheticians. The Manumed Comfort offers them an unrivalled examination couch in their practice or studio.

- Extra large back section makes this product the ideal examination couch
- Easy change from sitting to lying position with gas spring support
- Wheel raising mechanism for easy movement of the couch
- Height can be adjusted hydraulically (no mains lead needed) or electrically (very convenient, smooth and quick)
- Independently adjustable leg sections (optional)









*"Comfort and convenience make this Manumed the ideal examination couch."*

#### MANUMED COMFORT 3-SECTION




Type	
612	Manumed Comfort, hydraulic H/L
622	Manumed Comfort, electric H/L
632	Manumed Comfort with electrical section adjustment, electric H/L

#### MANUMED COMFORT 3-SECTION WITH SEPERATE LEG SECTIONS



Type	
613	Manumed Comfort with seperate leg sections, hydraulic H/L
623	Manumed Comfort with seperate leg sections, electric H/L


#### MANUMED COMFORT 2-SECTION




Type	
312	2-section, back-section, hydraulic H/L
322	2-section, back-section, electric H/L

#### TECHNICAL SPECIFICATIONS

Size tabletop (l x w):	200 x 67 cm
Height adjustment:	45-95 cm
Lifting time (min.-max.):	approximately 25 sec.
Lifting capacity:	150 kg
Force hydraulic pump:	10.000 N
Force electrical motor:	6.000 N
Power supply:	100/120/230V AC 50/60 Hz
Current consumption:	1.0 A max.



#### FRAME COLOURS



#### UPHOLSTERY COLOURS

012	013	004	169	141	199
155	138	019	139	127	135
811	812	852	888	933	851




## Manumed Exercise couch

Neurological exercise therapy with the Manumed Bobath. The large, extra wide treatment surface of 100 cm to 120 cm makes it easy to support and assist patients during therapeutic exercises. The Bobath and Vojta motion methods can be applied to perfection owing to the unmatched stability. A sturdy couch affording both patient and practitioner the most comfortable treatment positions.

- Very large surface (100 or 120x200 cm)
- Wheel raising mechanism for easy movement of the couch
- Electrical or hydraulic height adjustment
- The 2-section model can be used as a neurological exercise or treatment couch



### MANUMED EXERCISE 2-SECTION



#### Type


511	Manumed Exercise 2-section, hydraulic H/L 43 - 93 cm
521	Manumed Exercise 2-section, electric H/L 43 - 93 cm

### TECHNICAL SPECIFICATIONS


<b>Size tabletop (lxw):</b>	200 x 100 cm for models 500, 510 and 520 (optional 120 cm width for 510 and 520)
<b>Size tabletop (lxw):</b>	200 x 120 cm for models 511 and 521
<b>Height adjustment:</b>	43-93 cm
<b>Tijdsduur hoogte-</b>	
<b>Lifting time (min.-max.):</b>	approximately 25 sec.
<b>Lifting capacity:</b>	150 kg
<b>Force hydraulic pump:</b>	10.000 N
<b>Force electrical motor:</b>	6.000 N
<b>Power supply:</b>	100/120/230V AC 50/60 Hz
<b>Current consumption:</b>	1.0 A max.



### MANUMED EXERCISE 1-SECTION HEIGHT ADJUSTABLE




### MANUMED EXERCISE 1-SECTION FIXED HEIGHT



#### Type

500	Manumed Exercise 1-section, fixed height van 55 cm
510	Manumed Exercise 1-section, hydraulic H/L 43 - 93 cm
520	Manumed Exercise 1-section, electric H/L 43 - 93 cm

### FRAME COLOURS



### UPHOLSTERY COLOURS

012	013	004	169	141	199
155	138	019	139	127	135
811	812	852	888	933	851




# Manumed Traction

For the application of lumbar and cervical traction therapy there is the Manumed Traction couch. The roller-mounted surface prevents friction below the hips and legs, so that the full traction force set is actually applied in the vertebral treatment. The Manumed Traction is available in two variants: with a fixed or adjustable working height. The variable-height traction couch also doubles as an ideal treatment couch, because the rolling sections can be locked in place.

- Rolling surfaces (can be locked in place)
- With mountings for the traction frame
- Fixation rails for belts come standard
- For use as a traction as well as a regular treatment couch





**MANUMED TRACTION HEIGHT ADJUSTABLE**



Type	
916	Manumed Traction, hydraulic H/L
926	Manumed Traction, electric H/L

## MANUMED TRACTION FIXED HEIGHT



Type
004 Manumed Traction, fixed height

## TECHNICAL SPECIFICATIONS

<b>Size tabletop (lxw):</b>	203 x 67 cm for type 916 en 926
<b>Size tabletop (lxw):</b>	195 x 67 cm for type 004
<b>Height adjustment:</b>	45-95 cm for type 916 en 926
<b>Fixed height:</b>	83 cm for type 004
<b>Lifting time (min.-max.):</b>	approximately 25 sec.
<b>Lifting capacity:</b>	150 kg
<b>Force hydraulic pump:</b>	10.000 N
<b>Force electrical motor:</b>	6.000 N
<b>Power supply:</b>	100/120/230V AC 50/60 Hz
<b>Current consumption:</b>	1.0 A max.

## FRAME COLOURS

11/11/2019

## UPHOLSTERY COLOURS

012	013	004	169	141	199
155	138	019	139	127	135
811	812	852	888	933	851



# MANUMED

## Manumed Tilt



A specially designed tilting table for moving patients steplessly from horizontal to vertical, between 0 and 90°. The ideal couch for neuromotor treatments, among others. Children can also be treated safely and without problems.

The Tilt comes as a 1-section or as a 2-section couch.


- Tilts from 0 to 90°
- Suitable for children
- Electrical height adjustment
- Large wheels for safe movement



### MANUMED TILT 2-SECTION




### MANUMED TILT 1-SECTION



### TECHNICAL SPECIFICATIONS


<b>Size tabletop (lxw):</b>	202 x 67 cm
<b>Height adjustment:</b>	55 - 108 cm
<b>Lifting capacity:</b>	150 kg
<b>Force electrical motor:</b>	6.000 N
<b>Power supply:</b>	100 / 120 / 230V AC 50/60 Hz
<b>Current consumption:</b>	3.0 A max.



### Type

722	Manumed Tilt 1-section, electric H/L and electric tilting
723	Manumed Tilt 1-section with leg fixation rails, H/L and electric tilting
724	Manumed Tilt 2-section, electric H/L and electric tilting
725	Manumed Tilt 2-section with leg fixation rails, electric H/L and electric tilting

### FRAME COLOUR



### UPHOLSTERY COLOURS

012	013	004	169	141	199
155	138	019	139	127	135
811	812	852	888	933	851


# Manumed Basic

The basic high-quality couch for examination and treatment, with a fixed working height of 80 cm and simple adjustment of sections. Available in two variants, with the 2-section head section design having a negative angle to facilitate massage and manipulation.

- Fixed working height of 80 cm
- Extremely stable



MANUMED BASIC HEAD-SECTION



MANUMED BASIC BACK-SECTION



## TECHNICAL SPECIFICATIONS

**Size tabletop (lxw):** 195x67 cm  
**Fixed height:** 83 cm  
**Weight capacity:** 135 kg  
**Upholstery:** 18 different colours available  
**Colour code frame:** grey = RAL 7035, legs chrome

## UPHOLSTERY COLOURS

012	013	004	169	141	199
155	138	019	139	127	135
811	812	852	888	933	851


Type	
001	Manumed Basic 2-section, head-section with nose hole, fixed height
002	Manumed Basic 2-section, back-section, fixed height

# Portable Massage Couch

A practical and safe alternative for applying a treatment, wherever you are. Easy to transport and quick to set up.

- Weighs only 15 kg
- Load capacity up to 135 kg
- Attachable head section with nose hole
- 56 cm wide sections
- Aluminium frame with foldable legs
- Fixed height of 73 cm

PORTABLE MASSAGE COUCH



## TECHNICAL SPECIFICATIONS

**Size tabletop (lxw):** 210x56 cm  
**Fixed height:** 73 cm  
**Weight capacity:** 135 kg  
**Upholstery:** grey  
**Frame:** Aluminium  
**Carrying weight:** 15 kg

3450007 Portable couch with head section



## Colours


A matching Enraf-Nonius couch for every interior. As in other practices, a pleasant ambience is important in physiotherapy practices. Clients must feel comfortable in the room where they are being treated. As each practice is furnished differently and has its own colour scheme, we regard it as important for the treatment couches to fit in with your look and feel. That is why we can manufacture our treatment couches in several frame and upholstery colours. The imitation leather can be supplied with a rugged or smooth finish.

Our imitation leather upholstery contains no phthalates and is flame retardant – a unique combination that makes this material suitable for medical use!



## Options Examination and Treatment Couches<sup>(1)</sup>

	Manumed Optimal 3-section	Manumed Optimal 2-section	ManuXselect	Manumed Osteo	Manumed Comfort	Manumed Exercise	Manumed Traction	Manumed Tilt	Manumed Basic
Number of available upholstery colours	18	18	6	18	18	18	18	18	18
Number of available frame colours	3	3	1	3	3	2 <sup>(2)</sup>	2 <sup>(2)</sup>	1	1
Upholstery extra hard or extra thick	■	■	✗	■	■	✗	■	✗	■
Section width options	57 / 67 / 80 cm	57 / 67 / 80 cm	67 / 80 <sup>(3)</sup> cm	57 / 67 / 80 cm <sup>(4)</sup>	67 / 80 cm	100 / 120 cm <sup>(5)</sup>	67 cm	67 cm	67 cm
Ergo-Wings <sup>(6)</sup>	■	■	✗	✗	■ <sup>(7)</sup>	✗	✗	✗	✗
Nose hole long / triangular / none	S / ■ / ■	S / ■ / ■	✗ / S / ✗	S / ■ / ■	S / ■ / ■	S / ✗ / ✗	S / ■ / ■	✗ / ✗ / ✗ <sup>(8)</sup>	S / ✗ / ✗
Operation									
All-around operation	■	■	■	■ <sup>(9)</sup>	■	■	■	✗	✗
Footswitch	■	■	■ <sup>(A)</sup>	■	S	■	■	✗	✗
Handswitch	✗	✗	■ <sup>(B)</sup>	■	■	✗	■	S	✗
Tableswitch	✗	✗	✗	■	■	■	■	✗	✗
Wheels									
Wheel raising set	S	S	■	S	S	S	S	S	✗
Motor									
Extra heavy motor	■	■	■	✗	✗	✗	✗	✗	✗



■ = Optional possibilities  
✗ = Not possible on this model  
S = Standard on this model


- (1) = Note! Please order these options together with the couch because of modification of the couch during production  
(2) = Not possible with type 500 and not with type 004  
(3) = 80 cm model only available in colours 004, 013 and 199  
(4) = Type dependent (see Manumed Osteo on page 12)  
(5) = Type dependent (see Manumed Exercise on page 15)  
(6) = Ergo-wings not possible in combination with fixation rails  
(7) = Only with type 312 and 322  
(8) = Nose hole only with type 724, 725  
(9) = Standard with type 446, 448 and 449



## Stools

3xxx010	Stool, height adjustable 40-53 cm <sup>(1)</sup>
3xxx011	Stool, mobile, height adjustable 59-85 cm <sup>(1)</sup>
3xxx012	Stool, mobile, provided with footrest, height adjustable 59-85 cm <sup>(1)</sup>
3xxx013	Stool, mobile, anatomic seat, height adjustable 43-56 cm <sup>(1)</sup>
3xxx020	Stool, mobile, saddle seat, height adjustable 59-85 cm <sup>(1)</sup>
3496561	Set of 5 soft wheels for mobile stools, perfectly suited for smooth surface

(1) = Available in 18 different colours




## Rolls and cushions

3xxx001	Halfroll, flat bottom, especially for the knee <sup>(1)</sup>	27x16x18 cm
3xxx002	Halfroll, flat bottom, small, especially for the ankle <sup>(1)</sup>	28x9x11 cm
3xxx003	Halfroll, flat bottom, large <sup>(1)</sup>	56x16x18 cm
3xxx021	Halfroll, flat bottom <sup>(1)</sup>	50x18x10 cm
3xxx022	Halfroll, flat bottom <sup>(1)</sup>	56x25x15 cm
3xxx023	Massage roll <sup>(1)</sup>	50x Ø15 cm
3xxx004	Neck cushion, hourglass type <sup>(1)</sup>	30x15 cm
3xxx005	Massage cushion <sup>(1)</sup>	34x27x8 cm
3xxx006	Face cushion, with breathing hole <sup>(1)</sup>	40x20x15 cm
3xxx015	Wedge shaped cushion <sup>(1)</sup>	50x30x15 cm
3xxx016	Wedge shaped cushion <sup>(1)</sup>	50x40x30 cm
3xxx024	Wedge shaped cushion <sup>(1)</sup>	50x40x12 cm
3xxx017	Exercise roll <sup>(1)</sup>	100xØ30 cm
3xxx029	Plug for nosehole, small <sup>(1)</sup>	13x9 cm
1550813	Plug for nosehole, triangular shaped (only available in black)	16,5 x 21 cm
3445297	Wedge shaped, Kaltenborn, rubber model (only available in black)	22,5 x 10,5 x 8,5
3446161	Support bag with sand 1,2 kg <sup>(2)</sup>	15x25 cm
3446162	Support bag with sand 2 kg <sup>(2)</sup>	18x27 cm
3446166	Support bag with sand 4 kg <sup>(2)</sup>	25x35 cm
3446167	Support bag with sand 7 kg <sup>(2)</sup>	30x45 cm

(1) = Available in 18 different colours

(2) = Available in 1 colour only (grey 013)




# ENRAF-NONIUS


## Paper roll holders and paper rolls

		Manumed Optimal 3-section	Manumed Optimal 2-section	ManuXelect	Manumed Osteo	Manumed Comfort	Manumed Bobath	Manumed Traction	Manumed Tilt	Manumed Basic
3450104	Paper roll holder for Manumed Optimal, for 40, 50 & 60 cm rolls, mounted under foot-section	■	■							
3450105	Paper roll holder for Manumed Optimal, for 40, 50 & 60 cm rolls, mounted under headsection	■	■							
3450066	Paper roll holder, for 40, 50 & 60 cm rolls, mounted on the lower frame				■	■		■		
3450065	Paper roll holder leg section for 60 cm rolls					■				
3450052	Paper roll holder back section, for 40, 50 cm rolls					■			■ (1)	
1550800	Paper roll holder ManuXelect for 40 & 50 cm rolls			■						
3445280	Paper roll stand (floor type) for 40, 50 & 60 cm rolls	■	■	■	■	■	■	■	■	■
3445284	Paper roll, width 40 cm, roll-length 150 m, per 5 rolls	■	■	■	■	■		■	■	■
3445278	Paper roll, width 50 cm, roll-length 150 m, per 5 rolls	■	■	■	■	■		■	■	■
3445279	Paper roll, width 60 cm, roll-length 150 m, per 4 rolls	■	■		■	■	■	■	■	■

(1) = See page 24



3450104



3450105



3450066



3450065



3450052



1550800



3445280



3445284



3445278




3445279

## Miscellaneous


		Manumed Optimal 3-section	Manumed Optimal 2-section	ManuXelect	Manumed Osteo	Manumed Comfort	Manumed Exercise	Manumed Traction	Manumed Tilt	Manumed Basic
3446261	Set of 3 exercise stools, 50 x 38 x 32 cm, 40 x 32 x 21 cm and 30 x 22 x 14 cm	■	■	■	■	■	■	■	■	■
VK012668	Protective mat (transparent) 60 x 40 cm	■	■	■	■	■	■	■	■	■
3013032	Arm rests (set of two) for Manumed Osteo				■					
1550805	Gynaecological set <sup>(1)</sup>	■	■	■						
1550810	Gynaecological Comfort set <sup>(1)</sup>	■	■	■						
3450069	Raising set (30 mm)				■	■	■	■		
xxxxxx	Safety switch <sup>(2)</sup>	■	■	■	■	■	■	■	■	■



3446261




VK012668




3013032



1550805



1550810



3450069




xxxxxx

(1) = For use of this set the couch should be equipped with fixation rails (or rod), see page 23 for more information.

(2) = The safety switch differs per different couch types.  
1550801 = for the ManuXelect with 250 kg motor  
1550803 = for the ManuXelect with 150 kg motor  
VK010782 = for the Manumed (all types)

## Couch covers

3445240	Couch cover for Manumed Collection, with a cutaway for the nose	blue
3445241	Couch cover for Manumed Collection, with a cutaway for the nose	bordeaux
3445242	Couch cover for Manumed Collection, with a cutaway for the nose	white
3445243	Couch cover for Manumed Collection, with a cutaway for the nose	lavender
3445244	Couch cover for Manumed Collection, with a cutaway for the nose	jade
3445245	Couch cover for Manumed Collection, with a cutaway for the nose	beige
3445246	Couch cover for the Manumed Optimal	royal blue
3445247	Couch cover for the Manumed Optimal	off-white
3445034	Couch cover for the Manumed Optimal	bordeaux
3445261	Couch cover for the Manumed Optimal	blue
3445262	Couch cover for the Manumed Optimal	jade
3445267	Couch cover for the Manumed Optimal	white
3445259	Elastic straps for couch cover with clip fittings at both ends, set of 3	



3445240 / 3445261

3445241 / 3445034

3445242 / 3445267

3445243


3445244 / 3445262

3445245

3445246

3445247

3445259








# ENRAF-NONIUS

## Fixation accessories

			Manumed Optimal 3-section	Manumed Optimal 2-section	ManuXelect	Manumed Osteo	Manumed Comfort	Manumed Exercise	Manumed Traction	Manumed Tilt	Manumed Basic
A	3450134	Fixation rails Manumed Optimal <sup>(3)</sup>	■	■							
B	3450113	Clamp for armpit/shoulder fixation horns and for gynaecological set <sup>(1)</sup>	■	■	■						
C	1550815	Fixation horn, black <sup>(2)</sup> (set of 2)	■	■	■						
D	3xxx031	Pelvic support, incl. fixation board <sup>(3)(6)</sup>	■	■							
E	1550806	Fixation rail for fixation horn shoulder, armpit and for gynaecological set			■ <sup>(9)</sup>						
F	1550807	Fixation rail for fixation horn armpit and gynaecological set			■						
G	1550808	Fixation rail for fixation horn shoulder			■ <sup>(9)</sup>						
H	3450055	Fixation rails for Manumed Comfort type 312/322 (set of 4) <sup>(3)</sup>					■ <sup>(4)</sup>				
I	3450056	Fixation rails for Manumed type 400 & 600 (set of 4) <sup>(3)</sup>				■ <sup>(8)</sup>	■ <sup>(5)</sup>				
J	3445007	Clamp for armpit/shoulder fixation horn <sup>(1)</sup>				■	■ <sup>(5)</sup>				
K	3445006	Clamp for fitting straps (2 required) <sup>(1)</sup>				■	■ <sup>(5)</sup>				
L	3xxx007	Shoulder and armpit fixation horn <sup>(6)</sup>				■	■ <sup>(5)</sup>				
M	3450068	Clamp chrome, for shoulder fixation horn 3xxx030(1). 2 clamps for 2 horns.				■	■ <sup>(5)</sup>				
N	3xxx030	Shoulder fixation horn for Manumed Osteo and Traction <sup>(6)</sup>				■			■		
O	3xxx019	Pelvic fixation for Osteo and Traction <sup>(3)(6)</sup>				■			■		



(1) = To use these clamps the couch should be equipped with fixation rails

(2) = Fixation support requires horn 3450113 for the Manumed Optimal or fixation horn 1550806, 1550807 or 1550808 for the ManuXelect

(3) = ATTENTION! Needs to be ordered together with the couch because of adjustments of the section during production

(4) = For Manumed Comfort 2-section only

(5) = For Manumed Comfort 3-section only

(6) = Available in different colours (xxx is the colour code, see page 19)

(7) = Not in combination with ergo wings

(8) = Not with type 446, 448 and 449

(9) = Not with ManuXelect 2-section back

## Fixation straps

3445023	Fixation strap with velcro fastener 120 cm
3445014	Fixation strap with velcro fastener 160 cm
3445015	Fixation strap with velcro fastener 240 cm
3445475	Fixation strap with safety lock 250 cm
3445474	Fixation strap with safety lock 350 cm



3445023 3445014 3445015 3445475 3445474

## Accessories Manumed Traction<sup>(1)</sup>

		Manumed Traction Fixed Height	Manumed Traction Height adjustable
3445021	Thoracic belt for Manumed Traction	■	■
3445022	Pelvic belt for Manumed Traction	■	■
3445020	Traction frame 480 mm width, for type 004	■	
3445038	Traction frame 430 mm width, for type 916/926		■
3445512	Mobile traction frame for cervical traction	■	
3xxx014	Flexistool with separate height adjustment <sup>(2)</sup>	■	■



3445021 3445022 3445020 3445038 3445512 3xxx014



## Accessories Manumed Tilt

		Manumed Tilt with central fixation rail, 1- or 2-section	Manumed Tilt without central fix. rail, 1- or 2-section
3496629	Battery pack	■	■
3445043	Arm rests (set of two)	■	■
3445044	Work table	■	■
3445065	Fixing bar for arm rests and/or the work table	■	■
3445015	Strap with back-to-back fastening, 240 cm	■	■
3445024	Fixation straps for the knees or thorax fixation	■	■
3445046	Fixation straps for the knees for the Manumed Tilt with central fixation rails		■
3450052	Paper roll holder for 2-section	■	■
3450070	Paper roll holder for 1-section	■	■
3451055	Set of wheels with 2-sided operation <sup>(3)</sup>	■	■



3496629 3445043 3445044 3445065 3445015 3445024



3445046 3450052 3450070 3451055



(1) = For additional traction accessories we refer to the documentation of the Eltrac 471

(2) = Available in different colours (xxx is the colour code, see page 19)

(3) = ATTENTION! Needs to be ordered together with the couch because of adjustments of the product during production

# WWW.ENRAF-NONIUS.COM



## PASSION FOR MOTION

Your distributor:

### ENRAF-NONIUS INTERNATIONAL

Vareseweg 127

3047 AT Rotterdam

P.O. Box 12080

3004 GB Rotterdam

The Netherlands

Tel.: +31-(0)10 - 203 06 00

Fax: +31-(0)10 - 203 06 99

E-mail: [info@enraf-nonius.nl](mailto:info@enraf-nonius.nl)

Website: [www.enraf-nonius.com](http://www.enraf-nonius.com)